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
## General Information

**Intended use.** This irrigation pump has been designed for pumping irrigation fluid into the abdomen during laparoscopic procedures.

Do not use this product for any other purpose.

Do not use for procedures in which a defined pressure value must be preselected and reached, e.g. in hysteroscopy or arthroscopy.

For more information regarding the application and safe use of this product, refer to the Olympus Endoscopy System Guide.

**Medical use.** This product is to be used only in a medical facility by trained medical personnel.

**Safety advisory.** Follow all warnings and cautions contained in this manual and in the Olympus Endoscopy System Guide.

**Review all instruction manuals thoroughly.** Before use, review this manual, the Olympus Endoscopy System Guide, and the manuals for all other equipment which will be used during the procedure.

### **WARNING!**

An insufficient understanding of the information in these manuals can result in serious injury or equipment damage.

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CE 0124

**Use of this manual.** This manual contains valuable specification, care, and problem solving information which will help you use the product safely and effectively.

This manual can be a useful training aid for new users of endoscopic equipment and is a good reference for experienced users. Keep it in a safe, accessible location.



**Medical Device Directive (Europe).**

This product complies with the requirements of Directive 93/42/EEC concerning medical devices.

**Upon Delivery**


**Inspect immediately.**

- Make sure that the shipping case contains all of the items listed below.
- Contact your Olympus representative or an authorized service center if any items are missing or damaged.

**Contents A5889.**

- Pump "Eco-Pump"
- Power cable A5130, German plug
- Mount A5888
- Instruction manual
- 2 spare fuses

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**Unpacking.** Remove the instrument from the packaging. Do not use the packaging for routine storage.

---

**Manufacturer/Distributor**

**Manufacturer of the Eco-Pump:**  
W.O.M. WORLD OF MEDICINE AG  
Alte Poststrasse 11  
D-96337 Ludwigsstadt

**Distributor of the Eco-Pump:**  
Olympus Winter & Ibe GmbH  
Kuehnstr. 61  
D-22045 Hamburg

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

## Terminology

### Front.

- ① Pictograph for insertion of tubing set
- ② Power switch
- ③ Start/stop key
- ④ Roller wheel
- ⑤ Pump head with pressure sensor
- ⑥ Switch for pressure chamber

### Back.

- ① Identification plate
- ② Mount adapter
- ③ Connector for potential equalization
- ④ Fuse holder
- ⑤ Connector for unit plug

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Service

## Symbols and Labels

Refer to the accompanying instruction manual.

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
## Before Use

### WARNING!

**Infection control risk.** Properly reprocess all equipment before first and each subsequent use following the instructions in this manual and in the Olympus Endoscopy System Guide. Improper and/or incomplete reprocessing can cause infection of the patient and/or medical personnel.

### WARNING!

**Explosion hazard.** Install and operate this unit outside of the anaesthetic gas zone of risk.


**Ambient conditions.** Install the unit on a flat, dry surface.

Temperature	10 °C - 40 °C
	50 - 104 °F
Relative humidity	30 % - 75 %.

## Mounting on a Stand

### NOTE!

**Pump position.** Positioning the pump on the same level as the patient is recommended. Make sure that the pump position does not make the stand unstable, i.e. the maximum mounting height for the pump should be e.g. 1.10 m if the stand's foot has a diameter of 42 cm.

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

### Mounting the pump.

- Attach the mount (A5888) to a suitable stand and secure it at the required height.
  - Slide the two pins of the mounting adapter on the pump's back into the two boreholes of the mount.
- The pump is now securely attached to the stand.


### Mains Connection

#### WARNING!

**Risk of electric shock.** Connect this unit only to a grounded AC mains supply that complies with IEC requirements.

Connect this unit to a supply circuit that corresponds with the mains voltage indicated in the fuse holder's window on the unit's back:

230 for the range of  $230\text{ V} \pm 10\%$

#### CAUTION!

**The unit is not splash-proof.** Protect the unit from fluid. Keep the power socket dry.

Should any fluid enter the unit or the socket, do not use the unit.

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### CAUTION!

**Risk of damage caused by incorrect voltage.** Make sure that the correct local voltage has been preselected before delivery. If not, contact an authorized service center. Otherwise, the unit will be damaged.



### Mains connection.

- Press the power switch on the unit's front to the "O" position (off).
- Plug the power cable into the back of the unit.
- Plug the other end of the cable into the mains socket.
- If local safety regulations require, connect the unit to a potential equalization line.

### Tubing Connection

#### Select a compatible tubing set.

- Tubing set A4054: Sterile, intended for single-use only. Do not use if the packaging is damaged. Do not attempt to reuse or reprocess.
- Tubing set A4055: Non-sterile, intended to be used several times. Sterilize before use.

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

**Terminology of the tubing sets.**

- ① Cannula
- ② Tubing clamp
- ③ Y-piece
- ④ Pump section
- ⑤ Pressure chamber
- ⑥ Male Luer-lock adapter


**CAUTION!**

**Risk of damage.** Make sure that the tubing set is connected without being kinked or pinched.


**CAUTION!**

**Risk of damage.** Do not insert the tubing set's pump section while filled. Otherwise the pressure chamber's membranes may be damaged.

**Insert the tubing set's pump section.**

- Before first use, remove the light-blue transport safeguard from the pressure chamber.
- Insert the pressure chamber into the lower receiver of the pump head. The pressure chamber must click into position if correctly inserted.
- Stretch the tubing set's flexible pump section around the roller wheel and hook onto the upper receiver of the pump head.

**WARNING!**

**Irrigation fluid.** The irrigation fluid must comply with national pharmaceutical regulations.

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### CAUTION!

**Risk of implosion.** Use only flexible irrigation bags.  
Do not use glass bottles.

### Connect irrigation bags.

- Position two irrigation bags containing an appropriate irrigation fluid on a suitable stand at a distance of one meter above pump level.
- Close the tubing clamps.
- Insert the cannulas into the irrigation bags.

### Connect a suction/irrigation tube.

- Select a compatible suction/irrigation tube (with valves), e.g. A5605 or A5796 + A5797/98/99.
- Suction/irrigation tube with Luer-lock connector: Connect the tubing set's Luer-lock connector to the suction/irrigation tube's irrigation port.
- Suction/irrigation tube without Luer-lock connector: Detach the Luer-lock connector from the tubing set and connect the tubing set directly to the suction/irrigation tube's irrigation port.
- Connect the suction/irrigation tube's suction port to the central operating room suction system.
- Open one or both tubing clamps.

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## Use

**Before use.** Make sure that the unit and accessories have been properly reprocessed, inspected, and tested.

### Functional Inspection



**Functional inspection.** Prior to each use a functional inspection as described below must be performed.

#### NOTE!

Irrigation can be started or stopped with the suction/irrigation tube's valve. The irrigation flow rate, however, cannot be controlled.

#### Switch the unit on.

- Press the power switch to "I" (on). The LED lights up. A signal tone sounds for 1 s and indicates an automatic self-test of the unit.

#### Test function.

- Close the suction/irrigation tube's valve.
- Press the start/stop key. The LED lights up. A prepressure of 400 mm Hg builds up in the tubing system in order to provide full irrigation capacity.
- Open the suction/irrigation tube's valve. Irrigation starts.
- Close the suction/irrigation tube's valve. Irrigation stops. The prepressure for full irrigation capacity builds up again.

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- Press the start/stop key again. The LED goes out. The pressure in the tubing system is relieved automatically (the roller wheel rotates backwards).
- Press the start/stop key once more. The LED lights up. The prepressure builds up again. Functional inspection is completed. The pump is ready for use.


#### **Result of functional inspection.**

Use the unit only if the result of the functional inspection is positive.

The inspection result is negative, if:

- the inspection is not performed as described.
- the LED is defective.

If this is the case, do not use the unit!  
Contact your Olympus representative or an authorized service center.

#### **Procedure**

##### **Displace air from the tubing system.**

- Place the suction/irrigation tube's distal end into a receptacle.
- Open the suction/irrigation tube's valve: Displace all air bubbles by letting irrigation fluid flow through the tubing system and the suction/irrigation tube into the receptacle prior to each use.

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
### Introduce the suction/irrigation tube.





- Introduce the suction/irrigation tube into the abdomen through a compatible trocar tube. Refer to the suction/irrigation tube's instruction manual.

### Replacement pump.

- Prepare a second pump and accessories for replacement in case of failure/malfunction of the pump and/or accessories.

### Intervention.

- Open the suction/irrigation tube's valve to start irrigation.
- Perform the intervention.

### After use

#### Finish the intervention.

- Withdraw the suction/irrigation tube from the trocar tube.
- Press the start/stop key. The LED goes out.
- Press the power switch to "O" (off).
- Remove the tubing set from the pump.

#### Single use tubing set A4054.

- Dispose of the tubing set after use. Do not attempt to reuse or reprocess!

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**Reprocessing.**

- Bring the instrument to the reprocessing area.
- Reprocess all equipment following the instructions in this manual and in the Olympus Endoscopy System Guide.

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## Repair

**Policy.** Contact your Olympus representative or an authorized service center for repair and warranty information.

### **DANGER!**

**High voltage.** Do not open the unit. No user-serviceable parts inside.

**Repair guidelines.** Follow the general repair guidelines in the Olympus Endoscopy System Guide.

## Disposal

### **Waste electrical and electronic equipment.**

In accordance with European Directive 2002/96/EC on Waste Electrical and Electronic Equipment, this symbol indicates that the product must not be disposed of as unsorted municipal waste, but should be collected separately.

Refer to your local Olympus distributor for return and/or collection systems available in your country.

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
## Two-Yearly Inspection

**Inspection date.** An authorized service center must inspect this unit every two years.

A sticker on the back of the unit specifies the next inspection date.

**Inspection certificate.** The owner of this unit must make sure to receive a certificate from the authorized service center for every maintenance and/or repair. This certificate must contain:

- type and scope of the service
- service date
- name of the service company
- signature

## Annual Inspection

**Annual inspection test.** A hospital technician must inspect this unit every twelve months. This annual inspection includes a safety test and a basic functional test.

### Test record.

- Record the results of the annual inspection.
- Copy the form at the end of this instruction manual.
- Keep the test records together with this instruction manual.

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## **Safety Test**

### **Visual inspection.**

Visually inspect the unit. Make sure that:

- the fuse corresponds with the specifications (refer to the "Technical Data" section),
- the labeling on the unit is legible,
- the mechanical condition of the unit does not impair safe use,
- no safety-relevant contamination of the unit exists.

### **Measurement of the leakage current.**

- Measure the leakage current in accordance with DIN EN 60601-1/IEC 601-1.

### **Measurement of the protective conductor resistance.**

- Connect the power cable.
- Measure the protective conductor resistance in accordance with DIN EN 60601-1/IEC 601-1.  
The maximum value is 0.2  $\Omega$ .

### **Measurement of the insulation resistance.**

- Measure the insulation resistance with 500 - 700 V DC.  
The minimum value is 50 M $\Omega$ .  
Measurement of electric strength using high voltage is not allowed.

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

## Basic Functional Test

### Preparation.

- Make sure that the power cable is connected.
- Insert the tubing set and connect two irrigation bags (refer to the instructions under "Tubing Connection").
- Do **not** connect a suction/irrigation tube.
- Place the end of the tubing set into a graduated receptacle (graduation up to 0.2 l).
- Press the power switch to "I" (on). The LED lights up.
- Press the start/stop key. The LED lights up. Irrigation starts.
- Make sure that the tubing system is completely filled with irrigation fluid and no air bubbles remain.
- Clamp the tubing set at the instrument outlet.
- Empty the receptacle.

### Perform the test.

- Open the tubing set at the instrument outlet and take the time for 1 minute.
- Press the start/stop key after 1 minute: The level of irrigation fluid in the receptacle must be 1.8 l ( $\pm 10\%$ ).

**Result.** The test result is positive if the value given above is reached.

In case the given value is not reached contact an authorized service center.

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## Troubleshooting

### **Start/stop key flashes, three short warning tones sound.**

*Problem:*

The pressure chamber is not correctly inserted into the pump head. Irrigation does not start.

*Countermeasure:*

- Remove the pressure chamber from the pump head.
- Reinsert the pressure chamber into the pump head. The pressure chamber is correctly inserted if it clicks into position audibly.
- If the problem cannot be solved, contact an authorized service center.

### **Roller wheel stops.**

*Problem:*

The actual pressure has exceeded 500 mmHg for more than 5 seconds.

*Countermeasure:*

- Reduce the intraabdominal pressure. Irrigation starts again when the intraabdominal pressure has fallen below 500 mmHg.
- If the problem cannot be solved, contact an authorized service center.

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### **Continuous warning tone sounds.**

#### *Problem.*

An error has occurred at the pressure sensor or in the pressure measuring electronics.

#### *Countermeasure:*

- Switch off the unit and, if necessary, remove the pressure chamber.
- Switch the unit on again after 10 seconds.
- If the problem cannot be solved, contact an authorized service center.

### **Fuse Replacement**


#### **Power switch does not illuminate.**

#### *Problem.*

The power switch does not illuminate after switching the unit on.

#### *Countermeasure:*

- Make sure that the unit has been correctly connected to the power supply.
- Examine the unit's fuses and replace if necessary.
- If the power switch still does not illuminate, contact an authorized service center.

#### **Remove fuse holder.**

- Remove the power cable from the unit.
- Remove the fuse holder next to the connector for the unit plug on the unit's back.
- Replace the fuses.  
The standard set contains two spare fuses.

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**Reinsert fuse holder.**

- Reinsert the fuse holder.
- Reconnect the power cable.
- Switch the unit on.

**Disposal**

**Disposal of the unit.** If the unit and/or parts of it are defective, dispose of them according to the relevant national regulations.

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