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
AMULET Innovality

*AMULET Innovality*

FUJIFILM DIGITAL MAMMOGRAPHY SYSTEM

FROM TECHNOLOGY INTO INNOVATION. FROM INNOVATION INTO CARE

# Lasting smiles for women worldwide

Innovation and quality in mammography


AMULET Innovality—the result of Fujifilm’s ongoing “innovation” and commitment to providing top “quality” mammography services. The Innovality utilises Fujifilm’s unique a-Se direct conversion flat panel detector (FPD)\* to produce clear images with a low X-ray dose. This system makes use of intelligent AEC (i-AEC) combined with a new image analysis technology to automatically optimize the X-ray dosage for each breast type. AMULET Innovality is a highly advanced mammography system which offers an extremely fast image interval of just 15 seconds. With this system, Fujifilm furthers the provision of high quality examinations with superior image diagnostic accuracy.


\*Using a HCP (Hexagonal Close Pattern) TFT array

Origin of the name: Origin of the name: With its mammography solutions Fujifilm hopes to be an “Amulet” — always there to protect women’s health and allow them to be true to themselves, vibrant and beautiful. The AMULET series aims to provide top-class digital mammography solutions that can be customised to meet every sites needs.


## Unique detector for fast, low dose examinations

AMULET Innovality employs a direct-conversion flat panel detector made of Amorphous Selenium (a-Se) which exhibits excellent conversion efficiency in the mammographic X-ray spectrum. The new HCP (Hexagonal Close Pattern) detector efficiently collects electrical signals converted from X-rays to realize both high resolution and low noise. This unique design makes it possible to realize a higher DQE (Detective Quantum Efficiency) than with the square pixel array of conventional TFT panels. With the information collected by the HCP detector, AMULET Innovality creates high definition images with a pixel size of 50  $\mu\text{m}$ , the finest available with a direct-conversion detector.

Conventional square pixel

AMULET Innovality hexagonal pixel

This low-noise and high-speed switching technology allows tomosynthesis exposures with a low X-ray dosage and short acquisition time to be performed. Fast image display is also possible, realizing a smooth mammography workflow from exposure to image display.

## ISC – Optimized contrast and low X-ray dose using a Tungsten Target

Image-based Spectrum Conversion\* (ISC) technology can be used to optimize contrast in an image. ISC analyzes images to compensate for variations in contrast due to the density of mammary glands, amount of fat and X-ray spectrum. ISC aims to ensure that images display adequate contrast even with the use of a high energy, low-dose X-ray beam. This technology allows sites that previously exploited the superior contrast of a Molybdenum target to realize the dose advantages offered by the use of Tungsten without having to compromise image contrast.

\*Based on image analysis the appearance is adjusted to emulate the image quality with the simulated "optimal" spectrum.


# High quality images for easier diagnosis

Tomosynthesis: making it possible to observe the internal structure of the breast

In breast tomosynthesis, the X-ray tube moves through an arc while acquiring a series of low-dose X-ray images. The images taken from different angles are reconstructed into a range of Tomosynthesis slices where the structure of interest is always in focus.

The reconstructed tomographic images make it easier to identify lesions which might be difficult to visualize in routine mammography because of the presence of overlapping breast structures.

The Tomosynthesis function on AMULET Innovality is suitable for a wide range of uses, offering two modes to cater for various clinical scenarios. Standard (ST) mode combines rapid exposure timing and efficient workflow with a low X-ray dose while High Resolution (HR) mode makes it possible to produce images with an even higher level of detail, allowing the region of interest to be brought into clearer focus.

## Breast screening and diagnostic mammography: two modes to answer different clinical needs

### Routine mammography and breast screening: ST mode

- Acquisition angle:  $\pm 7.5^\circ$
- Pixel size: 100/150  $\mu\text{m}$

The smaller angular range and fast image acquisition allow Tomosynthesis scans to be quickly performed with a relatively low X-ray dose.

### Diagnostic mammography: HR mode

- Acquisition angle:  $\pm 20^\circ$
- Pixel size: 100/50  $\mu\text{m}$

With a larger acquisition angle the depth resolution is improved. This allows the region of interest to be defined more clearly and brought into clearer focus.

## Two types of image processing pattern

Image Processing Pattern 1

Enhances spicula and calcifications while keeping maximum contrast for the viewing of masses within the glandular tissue.

Image Processing Pattern 2

Maximizes the visualization of fine spiculations and calcification

## ISR – Iterative Super-Resolution Reconstruction

### Excellent-m 3D

The tomosynthesis iterative super-resolution reconstruction (ISR) method has been specifically designed to support routine and screening application of ST-mode (15° angle) tomosynthesis acquisition.

#### 1. Reducing graininess of image in low-dose tomography


Each image is being processed to selectively suppress the patterns that do not exist in human body architectures (as the graininess related to noise), to reduce distractive details in the event of low-dose acquisition.

#### 2. Suppressing interference of human body architectures at different depths

In the process of reconstructing the breast structure from multiple 2D images, details as calcifications, masses, spiculas, mammary gland and other signals that emerge from different depths are removed to reproduce the correct structures at the focus depth with greater fidelity and accuracy.

#### 3. Restoring the fine-structure

Our super-resolution technology is introduced to restore the fine-structure of microcalcifications and other details, the visibility of which might be impaired by the movement of the X-ray tube.

## S-View function for synthesized 2D images reconstruction

S-View function allows reconstruction of synthesized 2D images from the tomosynthesis acquisition data: both 2D and 3D images are provided for a more comprehensive interpretation.


## Easy operation and patient comfort

### Intelligent AEC optimizes the X-ray dose for each breast type

Intelligent AEC has advantages in defining the optimal dose for an examination compared to conventional AEC systems where the sensor position is fixed.

Through the analysis of information obtained from low-dose pre-shot images, Intelligent AEC makes it possible to consider the mammary gland density (breast type) when defining the x-ray energy and level of dose required.

Able to be used even in the presence of implants; intelligent AEC enables more accurate calculation of exposure parameters than is possible with conventional AEC systems. By allowing the use of automatic exposure for the implanted breast, Intelligent AEC can further enhance examination workflow.


### Patient information display

The information shown on the display **A** at the base of the exposure unit can be switched between patient information (ID, name, date of birth, etc.) and positioning information (angle of swivel arm, compression force and breast thickness). Positioning information can also be confirmed on the display **B** on the compression arm.

### Dedicated mammography AWS (Acquisition Workstation)

#### Optimal examination workflow

- Integrated X-ray controller allows setting and confirmation of exposure conditions on a single screen.
- Examination screen can be split and switched between 1, 2, or 4 image display.
- Individual images can be immediately output to a PACS, viewer or printer during an examination.
- Density and contrast can be easily adjusted while viewing images.
- Alignment of left and right images can be adjusted both automatically and manually.

#### High definition second monitor (3M/5M: Optional)

- A second, high resolution monitor can be added to the AWS making it possible to display previous images recalled from a PACS to ensure the mammographer has access to previous images at all times.
- For Tomosynthesis, reconstructed images can be displayed.

#### Variable height operator desk

- A new desk has been designed to maximize operators comfort and system ergonomics, which is featured with an automated adjustable height mechanism: each user is able to operate in the best conditions and safety.

## AMULET Harmony

AMULET Harmony incorporates a range of mammography solutions specifically designed to maintain a harmonious examination environment and foster an atmosphere of trust between mammographers and their patients.

### Mood lighting to ease patient anxiety

Warm indirect lighting is used to illuminate the exposure stand, helping patients to relax and allowing examinations to be performed with minimal stress.

### Decorative labels adaptable to each room environment

Five different stand labels are available to add a gentle ambience. Each site can choose a stand appearance that best suits the examination environment, thus relieving patient stress and anxiety.

### Compression Paddles

Fujifilm Innovality features a wide lineup of compression paddles, in order to guarantee the best patient and breast positioning in all conditions. Special sizes and shapes are available as well as the default formats 18x24 and 24x30.

### Fit Sweet Paddles

This type of compression paddle (available in 18x24 and 24x30 format) fits to the shape of the breast, allowing pressure to be evenly applied while holding the breast securely and ensuring the breast tissue is adequately separated. Models with the lateral shift function (see below) are also available in the lineup.

### Shift Compression Paddles

This small compression paddle (18x24) can be positioned in the middle, right or left side of the detector at any time of examination according to the positioning of the patient.

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## EC Declaration of Conformity

13-CEM-02-U

Manufacturer: FUJIFILM Corporation  
Address: 26-30, Nishiazabu 2-chome, Minato-ku,  
Tokyo 106-8620, JAPAN  
Authorized Representative: FUJIFILM Europe GmbH  
Address: HEESENSTRASSE 31  
40549 DUESSELDORF, GERMANY  
Product(s): FUJIFILM Digital Mammography System  
AMULET Innovality  
Model No.: FDR MS-3500  
UMDNS Code: 20792  
GMDN Code: 37672  
Applicable Product Units: Serial 86621947 or later  
Classification(MDD, Annex IX): Class II b (Rule 10 )

We herewith declare in our own responsibility that the above mentioned product(s) meets the provisions of the following EC Council Directive(s) and standards.

Directive :  
Medical Device Directive : 93/42/EEC and their Annexes  
RoHS Directive : 2011/65/EU

Standards:  
MDD: EN ISO13485: 2012/AC:2012  
EN ISO14971:2012  
EN 60601-1:2006/AC:2010  
EN 60601-1-3:2008/AC:2010  
EN 62304:2006 /AC:2008  
EN 60601-1-6:2010  
EN 62366:2008  
EN 60601-2-45:2011  
EN ISO 10993-1:2009/AC:2010  
EN 60601-1-2:2007/AC:2010

RoHS: EN 50581:2012

EC Certificate for Directive 93/42/EEC: G1 17 07 20011 043  
Conformity Assessment Procedure: Annex II excluding (4)  
Notified Body: TÜV SÜD Product Service GmbH (Notified Body Number 0123)  
Ridlerstr. 65, 80339 München, Germany

Place and Date of issue

Kanaqawa, JAPAN

2018-4-11

Signature : 

Name : Naotake Mitsumori  
Function: Senior Engineering Scientist  
Quality Assurance and Regulatory Affairs Division  
Medical Systems Business Division  
FUJIFILM Corporation

## Annex I to EC Declaration of Conformity

Quality Assurance and Regulatory Affairs Division  
Medical Systems Business Division  
FUJIFILM Corporation

Please refer to the following additional information on EC Declaration of Conformity. This information is to be maintained with corresponding EC Declaration of Conformity.

EC Declaration of Conformity No. : 13-CEM-02-U

Material Description(s) applicable to  
EC Declaration of Conformity :

DR 3500 W 24X30 + E	NEEDLE HOLDER M 8G E
COMP PLATE 18X24 LOW HA A E	NEEDLE HOLDER EN E
COMP PLATE 18X24 HIGH HA A E	NEEDLE HOLDER H EV E
COMP PLATE 18X24 SMALL HA A E	NEEDLE GUIDE M/V E
COMP PLATE 18X24 FLEX HA E	NEEDLE GUIDE C 12G E
COMP PLATE 18X24 FS HA E	NEEDLE GUIDE C 14G E
COMP PLATE 18X24 S FS HA E	NEEDLE GUIDE C 16G E
COMP PLATE 18X24 SP HA E	NEEDLE GUIDE C 18G E
COMP PLATE 18X24 SP SMALL HA E	NEEDLE GUIDE C 20G E
COMP PLATE 18X24 FS SP HA E	NEEDLE GUIDE C 22G E
COMP PLATE 24X30 LOW HA A E	NEEDLE GUIDE EN E
COMP PLATE 24X30 SMALL HA A E	NEEDLE GUIDE TEST E
COMP PLATE 24X30 SSMALL HA E	NEEDLE GUIDE COVER E
COMP PLATE 24X30 FLEX HA E	NEEDLE GUIDE COVER 10 E
COMP PLATE 24X30 FS HA E	BIOPSY EQUIPMENT BOX E
COMP PLATE SPOT 9X9 HA E	LATERAL BIOPSY KIT HA E
COMP PLATE MAG1.8 16X20 HA E	LATERAL NEEDLE STOPPER M E
COMP PLATE MAGSPOT 9X9 HA E	LATERAL NEEDLE STOPPER C E
COMP PLATE ROUND SPOT HA E	LATERAL NEEDLE STOPPER V E
COMP PLATE ROUND MAG SPOT HA E	LATERAL NEEDLE STOPPER E E
COMP PLATE AXILLA 8X20 HA E	LATERAL NEEDLE STOPPER TEST E
COMP PLATE 2D BIOPSY HA E	ATTACHMENT FOR M RE E
COMP PLATE BIOPSY HA E	LATERAL ATTACHMENT FOR M RE E
COMP PLATE LATERAL BIOPSY HA E	DR 3500 TOMO KIT V2.0 E
DR 3500 ADD COLLI CONT KIT E	DR3500 TOMO S VIEW SOFT E
MAG TABLE HA E	DR3500 TOMO EXC SOFT E
EXPOSURE SWITCH REMOTE AR E	STEREO EXPOSURE KIT AR E
EXPOSURE SWITCH FOOT AR E	DR3500 S UP-GRADE KEY E
PROTECT BOARD FOR CL E	TOMO BIOPSY SOFT HA E
COVER FOR PROTECT BOARD E	FOOTSW SEQ MODIFY KIT A HA E
MONITOR ARM E	FOOTSW SEQ MODIFY KIT B HA E
DR 3500 TOMO FACE GUARD E	DR 3500 AC CORD TRANS 11M E
DR3500 GRIP UNDERTABLE KIT E	1SHOT PHANTOM 24X30 HA E
CHEST WALL PAD E	DR 2000 MBL PROTECTBD KIT AR E
AXILLA PAD LEFT E	DR 2000 MBL AWS PC FIX AR E
AXILLA PAD RIGHT E	DR 2000 MBL MONITOR FIX AR E
SLIP PROOF PLATE E	DR 3500 MBL AWS PC FIX E
PROTECT SHEET E	DR 3500 CEDM SOFT-L E
BIOPSY SPACER E	DR 3500 COLLIMATOR FOR CEDM E
POSITIONER KIT BIOPSY HA E	DR 3500 MBL KIT HA E
NEEDLE HOLDER M E	DR 3500 FACE GUARD T COMFORT E
NEEDLE HOLDER C E	DR 3500 CEDM MANUAL E
NEEDLE HOLDER V E	AWS PC COVER E
	DR 3500 COMFORT COMP SOFT-L E

See Annex to Declaration of Conformity of FDR-3000AWS (13-CEM-01-P or later)

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
## Delegation of Power of Attorney


### FUJIFILM Europe GmbH

the registered office at Heesenstraße 31, 405  
49 Düsseldorf, Germany  
corporate identification number (HRB): 10639,  
registered by the Regional Court in  
Düsseldorf,

operating in the Czech Republic through its  
branch:

**FUJIFILM Europe GmbH, organizační  
složka**, with its registered office in Prague 3,  
Žižkov, U Nákladového nádraží 1949/2, Postal  
Code: 130 00, the Czech Republic, ID number:  
24660736, registered in the Commercial  
Register maintained by the Municipal Court in  
Prague, Section A, Insert No. 71605,  
represented by Mr. Pawel Jerzy Augustyn,  
Manager of the Branch („Principal“).

Hereby grants a power of attorney to:

(„Representative“)

To represent the Principal in any legal and  
factual actions related to the Principal's  
business activities in so far as these may be  
necessary for proceedings pursuant to to  
Czech Act No. 134/2016 Coll., on Public  
Procurement, up to the total value of the  
tender of EUR 300,000 („Matter“).

Based on this power of attorney, the  
Representative is authorized to sign and  
execute on behalf of the Principal all  
documents and to perform further acts or  
negotiate as the Representative deems  
suitable or necessary for or in relation to the  
Matter.

Base on this power of attorney, the  
Representative is authorized, among others,  
to submit a bid, apply for a participation in a  
short list of bidders, take part in a tender  
dialogue, and any other acts necessary for the

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**FUJFILM Europe GmbH**
zastoupená organizační složkou FUJIFILM Europe GmbH, organizační složka