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
AGFA DR 600 - Equipment List	
1.01	3-Ray Table
1.02	RAM Table
1.03	Setting (Support)
1.04	Control Unit
1.05	Agfa Workstation
1.06	Control Unit and Workstation
1.07	Desk with chair and equipment
1.08	DR 600
* All dimensions without special	

### Unit arrangement and Site preparations

This was made in conformity with the architectural planning references and the project agreement. We reserve the right to make technical alterations. Before the construction work is completed it has to be ascertained that there are no changes in the original project.

This scope of delivery and services is described and determined in the order confirmation. All data serve for the on-site preparation for installation and setting of the DR 600 system in order to ensure correct operation.

This plan is not a construction drawing and is not to be used for carrying out construction work.

List of Documents	
Document No.	Description
01 FN Order - A107	3D System overview
02 FN Order - A102	Setup schedule
03 FN Order - A104	Site conditions - Floor
04 FN Order - A104	Site preparations - Ceiling
05 FN Order - A106	Site preparations - Side view
06 FN Order - A106	Site conditions - Workstation

**Room dimensioning**  
The indicated floor dimensions have to be observed on site. The geometry department has to be informed about possible deviations. Otherwise we cannot assume any guarantee for the successful implementation of the dimensions indicated in the planning documents.

**Dimensioning**  
All installation measurements apply to finished wall or floor and are to be checked prior to installing the unit.  
  
 ◆ Orientation point = reference point of the unit for planning and installation.


**Customer Release Acknowledgement**  
I agree to the preliminary site plan for implementation planning. We have been informed of the importance of meeting all rules and requirements. Agfa HealthCare assumes that all reader complies with on-site requirements and planning specifications, such as structural engineering, climate control, preparations for installation and removal, installation of the wall surface, the implementation and working of on-site preparations for installation of qualified specialists.

Customer	Project	Location	Date
FN Olomouc	AGFA	Olomouc	11.01.2011

**FN Olomouc**  
Olomouc  
Czech Republic

3D System overview	Room height finished floor to Concrete Ceiling
AGFA DR 600	Supporting Ceiling
	Scale: 1:25
	Unit: m
	1: _____
	Date: _____


Environment DR 600			
Temperature (indoor)	Operation	Transport	Storage
18°C to 21°C	18°C to 21°C	18°C to 21°C	18°C to 21°C
Relative humidity (max. condensation)	30% to 70%	30% to 70%	30% to 70%
Atmospheric pressure	970 hPa to 1013 hPa	970 hPa to 1013 hPa	970 hPa to 1013 hPa

Transport DR 600			
Box with labels	900x150x100 mm	1 x 20 kg	
Box with 19" Rack	900x250x100 mm	1 x 10 kg	
Box with Cooling Storage mode	900x250x120 mm	1 x 10 kg	
Box with Generator, 10k, 10k, 10k and 10k	900x250x100 mm	1 x 10 kg	
Box with Table, Table and 10k components (with 20kg Table)	900x250x100 mm	1 x 10 kg	

AGFA DR 600 - Equipment Legend	
1.01	1.01 1.01 Label
1.02	1.02 1.02
1.03	1.03 1.03
1.04	1.04 1.04
1.05	1.05 1.05
1.06	1.06 1.06
1.07	1.07 1.07
1.08	1.08 1.08
1.09	1.09 1.09
1.10	1.10 1.10

Measurements DR600	
Weight of transportable	Length of Transportable
Weight of transportable	Length of Transportable
Weight of transportable	Length of Transportable
Weight of transportable	Length of Transportable

**Electrical installation DR 600**

Generator Type	30 kW	65 kW	90 kW
Generator nominal rating	30 kW	65 kW	90 kW
max. P/N	50	65	80
max. I/A	300	330	400
max. V/E	150	150	100
Line voltage	400 VAC or 400V/430 VAC*		
Line frequency	50/60 ± 1Hz		
Max. line deviation	±10%, ±10%, @ 400V		
Line connection	3P+N+PE		
nominal capacitor wattage	24 kVA	28 kVA	44 kVA
Peak current sensitivity	30mA		
Fuse rating (slow blow) 400 input voltage	50A	50A	50A
Input current at max. output power rating per phase for 0.2 seconds	110A@400V 97A@430V	144A@400V 124A@430V	180A@400V 154A@430V
Max. line resistance at 400VAC per phase	0.30	0.20	0.120
Power factor Cos φ	0.8		
Max. conductive cross section of main lead-in at Transfer Point	25mm² for flexible wire 35mm² for solid wire		

- \* Temperature is available in two versions
- \* 400 V input
- \* With integrated pre-transformer which can be converted for 400V or 480V input

**Insufficient safety distances**

In accordance with the standard EN 61439, actions must be taken to ensure safety in danger zones A (danger zone) or in areas where the possibility of injury exists because of the position or design of stationary or moving parts of a technical unit.

These danger zones can be made safe through the use of protective facilities (barriers, covers, safety barriers such as light barriers, switching stops, switching mats, conductive warning signs hanging by wires).

Here, it must be ensured that the protective facilities are designed and positioned in such a way that the safety distances are not changed.

The safety distances at crushing points must, for example, be > 30 cm for the body, > 18 cm for the legs > 10 cm for the hands and > 12 cm for the feet and arms.

Notice: When setting up pieces of equipment and making other installations, which are not included in this planning, it must be ensured that the safety distances are maintained.

**Room lighting**

Ambient light in rooms where diagnostic tests place on image display devices (monitors) must meet the following requirements:

- Free of glare, controllable, reproducible setting of the lighting intensity (e.g. dimmer with scale)
- No reflections from windows, lamps and viewing boxes in the usual operating position of the image display devices.

This is a specification of DIN 6885 57 in Germany, which should also be complied with at other countries in regard to the lighting of rooms for diagnosis, imaging and treatment procedures. The intensity of the lighting is dependent on the type of procedure.

If only X-ray exposure are produced, the requirements for lighting for diagnostic imaging with image intensifiers apply (50 lx).

When images are displayed on monitors, the possibility of reducing the general lighting intensity must be provided (50 lx, if necessary, down to 1 lx).

Reflections and glare on the screen must be avoided (DIN 6973-3, EN 12464).

As a rule, therapy rooms require a general lighting with a nominal lighting intensity of 300 lx.

This also applies to rooms where patients are treated with physical, biological or electromagnetic procedures.

**Statics DR 600**

The system has to be installed on a solid surface with sufficient load carrying capacity. Wall or ceiling structure if the ceiling panels, e.g. suspended ceiling, do not have a sufficient bearing load, it must be removed and replaced by a concrete reinforcement slab (C20/25). If an appropriate system is not available, the ceiling must also be reinforced in the installation floor.

**Planning recommendations**

According to the German standard DIN 61813, an inter-connector and a device for measurement is required between the patient and the operator.

Consider the ceiling height, distance and light beams by arranging the lamps.

Darkening of windows in X-ray and control rooms is recommended.

**Notes on Preparations for installation**

The preparations for installation include routing of power lines as well as provision of suitable wiring and cables for transfer of digital information, setting up the inter-connector, installing floor boxes, switches in a main ground fault and emergency stop (off) switches, distribution boxes, transformers with separate windings and isolation monitors, installation of cable ducts, installation boxes and conduits according to plans completed in the same manner for installation and fastening of special low-voltage for mounting equipment on ceilings, independent lighting, floors and walls as well as installation of heavy load structures.

Complete for performing and supervising on-site installation preparation should be conducted with technically competent companies directly by the ordering party. The ordering party is responsible for safety and proper completion and supervision of all preparations for installation at the construction site in observance of all applicable legal regulations (e.g. X-ray regulations, radiation protection regulations) and all applicable general recognized rules of technology (e.g. VDE regulations, DIN standards).

Execution and supervision of installation preparations at the construction site and test observations of the standard operating conditions are not included in our offer.

The ordering party is responsible for checking the static characteristics and, where applicable, the air conditioning in the building to be equipped.

**Unit arrangement and Site preparations**

The work made in accordance with the construction planning references and the provided equipment. We reserve the right to make technical alterations. Before the construction work is completed it has to be ascertained that there are no changes in the original design.

The design of delivery and services described and determined in the order confirmation. All data serve for the general preparation for installation and setting of the DR 600 system in order to realize the final application.

The supply is not a production drawing and it is not to be used for carrying out construction work.

**Room dimensioning**

The indicated room dimensions have to be checked on site. The planning department has to be informed about possible deviations. Otherwise we cannot assume any guarantee for the accurate implementation of the dimensions indicated in the planning documents.

**Dimensioning**

If condition measurements apply to finished wall or floor and are to be checked prior to assembly, the unit.

Dimensions: 1000 mm, 1000 mm, 1000 mm

Dimensions: 1000 mm, 1000 mm, 1000 mm

Dimensions: 1000 mm, 1000 mm, 1000 mm

**Customer Release Acknowledgment**

I agree to the release of the data for information purposes. The data have been entered in the database of the AGFA DR 600 system and are available for all AGFA DR 600 systems. AGFA Healthcare does not check and make comparisons with other requirements and planning details given, such as national engineering or other control requirements for installation and design of installations, the will contact the engineering and installation of the site preparation for a solution in qualified companies.

Customer	Name	Address	City
AGFA	AGFA	AGFA	AGFA

**AGFA**

**FN Olomouc**  
**Olomouc**  
**Czech Republic**

**Setup Overview**

AGFA DR 600


Scale: 1:25

AGFA DR 600


# Site preparations Technical details - Side views


AGFA DR 600 - Equipment Legend	
1.01	X-Ray Table
1.02	Wall Stand
1.03	Casting Suspension
1.04	Generator off-line
1.05	Agfa Workstation
1.06	Generator on/off and exposure switch
1.07	Electromagnetic shielding
1.08	Di.D.F.F.F.

### Detailed sideview A

### Detailed sideview B

## Site Preparations Technical Details - Side views

We reserve the right to make technical alterations.  
 Details serve for the on-site preparation for installation and utility of the Agfa DR600 system in order to ensure correct operation.  
 This plan is not a construction drawing and is not to be used for carrying out construction work.

**Dimensioning**  
 All installation measurements apply to finished walls or floors and are to be checked prior to installing the unit.  
 Dimension lines + reference point of the unit for planning and installation.

**Room dimensioning**  
 The indicated room dimensions have to be checked on site. The planning department has to be notified about possible deviations. Otherwise we cannot assume any guarantee for the accurate implementation of the dimensions indicated in the planning documents.


**Customer Release Acknowledgement**  
 I agree to the performance on the basis of information provided. We have been informed of the importance of providing all views and requirements. Agfa Healthcare does not check and/or accept any compliance with site requirements and planning specifications, such as structural engineering, other building preparations for installation and electrical installations. We will collect the implementation and installation of site preparations for installation in specified areas only.

Customer	
First installation	
Customer approval	

**AGFA**  
 FN Olomouc  
 Olomouc  
 Czech Republic

Site preparations - Side views		Room height measured from to Concrete Ceiling Suspended Ceiling
AGFA DR 600	Scale: 1:25	
	Date:	
	Project:	

# Site preparations Technical details - Workstation


AGFA DR 600 - Equipment Legend	
1.01	Artex Table
1.02	Shel. Stand
1.03	Choking Support table
1.04	Utilization table
1.05	High. Workstation
1.06	Control unit and exposure switch
1.07	Used for work preparation
1.08	DR O.F.F.S
1.09	DR O.F.F.S


- 16x250-340V
- Network socket RJ45
- 600mm table with cable holder
- Deflection with mobile light
- Emergency stop with lock
- Dual start with outspooling
- Generator unit and exposure switch

1. Six sockets 16x250/340V, 500V, 16A
2. Network socket RJ45
3. Cable tray 150x30mm
4. Hole in the table with table holder
5. Generator unit and exposure switch
6. Dual start switch/outspooling
7. Emergency stop with lock
8. On-Off button with indicator light

## Detailed sideview C

## Detailed sideview D

## Site Preparations Technical Details - Workstation

We reserve the right to make technical alterations.  
 We do not assume any liability for the accuracy of the information provided in this document.  
 This document is a technical drawing and is not to be used for construction work.

**Dimensioning**  
 All dimension measurements apply to finished work or those who are to be checked prior to assembly.  
 Characteristic points = reference points of the unit for planning and installation.

**Room dimensioning**  
 The indicated room dimensions have to be checked on site. The planning department has to be informed about possible deviations. Otherwise we cannot assume any guarantee for the accurate implementation of the dimensions indicated in the planning documents.

**Customer Release Acknowledgement**  
 I agree to the preliminary as the basis for implementation planning.  
 We have been informed of the importance of handing all notes and requirements.  
 AGFA HealthCare does not check and monitor compliance with site requirements and planning specifications, such as structural engineering, climate control, preparations for installation and structural installations. We will contract the implementation and monitoring of on-site preparations for installation to a third party.

Customer	Name	Signature	Date
AGFA			

FN Olomouc  
 Olomouc  
 Czech Republic


Site preparations - Workstation

AGFA DR 600


Room Height (finished floor to concrete ceiling)	
Height	1.25
Scale	1:25

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