

# User-friendly Design

## Easy to Use Screen

## Easy Access to Patients Eyelids


The radical cut design allows direct access to patient eyelids.

## 3-D Auto Tracking and Auto Shot


The 3-D auto tracking and auto shot is faster, simpler, and more accurate than manual measurements.

Once alignment is completed, the measurement starts automatically.

## Joystick for Flexible Alignment


The joystick helps the operator make fine adjustments during alignment to improve the precision, even for eyes with poor fixation which cannot be tracked with automated tracking systems.


## Wireless LAN (WLAN)

The TONOREF™ III connects with PC and peripheral devices using wireless LAN (WLAN)\*, LAN cable, RS-232C cable, EyeCare Card, barcode scanner or magnetic card reader.

\*Available for products shipped for USA, Canada, and other countries that implement the R&TTE Directive.

## Automatic Anti Dew Heater

Automatic anti dew heater for measuring windows prevents condensation to provide accurate measurements in cooler rooms.

# TONOREF™ III Specifications

<b>Auto refractometer</b>	
Measurement range	Sphere -30.00 to +25.00 D (VD = 12 mm) (0.01 / 0.12 / 0.25 D increments) Cylinder 0 to ±12.00 D (0.01 / 0.12 / 0.25 D increments) Axis 0 to 180° (1° / 5° increments)
Minimum measurable pupil diameter	ø2 mm
Measurement area	ø1 to 6 mm
Chart	Scenery chart
<b>Auto keratometer</b>	
Measurement range	Curvature radius 5.00 to 13.00 mm (0.01 mm increments) Refractive power 25.96 to 67.50 D (n = 1.3375) (0.01 / 0.12 / 0.25 D increments) Cylindrical power 0 to ±12.00 D (0.01 / 0.12 / 0.25 D increments) Axis 0 to 180° (1° / 5° increments)
Measurement area	ø3.3 mm (R=7.7 mm), ø2.4 mm (R=7.8 mm)
<b>Non contact tonometer</b>	
Measurement range	1 to 60 mmHg (1 mmHg increments)
Measurement range setting	APC40, APC60 (APC=Automatic Puff Control), 40, 60
Working distance	11 mm
Eye fixation	Inner fixation light
<b>Non contact pachymeter</b>	
Measurement range	300 to 800 µm (1 µm increments)
IOP correction by corneal thickness	Automatic calculation
Retroillumination image	Available
Accommodation measurement range	0 to 10.00 D (0.01 / 0.12 / 0.25 D increments)
PD measurement range	30 to 85 mm (1 mm increments) (Near point PD: 28 to 80 mm at WD = 40 cm)
Corneal size measurement range	10.0 to 14.0 mm (0.1 mm increments)
Pupil size measurement range	1.0 to 10.0 mm (0.1 mm increments)
Auto tracking	X-Y-Z directions
Auto shot	Available
Display	Tiltable 7.0-inch color LCD with touch panel
Printer	Thermal line printer with easy loading and auto cutter
Interface	RS-232C: 2 port, LAN: 1 port, USB: 1 port Wireless LAN (WLAN)*: 1ch
Power supply	AC 100 to 240 V, 50 / 60 Hz
Power consumption	100 VA
Dimensions / Mass	260 (W) × 495 (D) × 505 (H) mm / 22 kg at ARK standard mode 260 (W) × 495 (D) × 460 (H) mm / 22 kg at NT standard mode 10.2 (W) × 19.5 (D) × 19.9 (H) " / 48 lbs. at ARK standard mode 10.2 (W) × 19.5 (D) × 18.1 (H) " / 48 lbs. at NT standard mode

\*Limited to the USA, Canada, and other countries that implement the R&TTE Directive.

Product / Model name: AUTO REF / KERATO / TONO / PACHY METER TONOREF III

Specifications may vary depending on circumstances in each country.

Specifications and design are subject to change without notice.

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