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Biological Microscope

**CX23**

Bringing Reliability to the Classroom

# User-Friendly Design and Outstanding Optical Performance

Ergonomic grip

- Ergonomic grips add safety when retrieving the microscope from high places.
- Colored grips indicate the appropriate places to hold the microscope.
- Smooth, rounded design eliminates sharp edges.
- The CX23 microscope is one of the lightest in its class, with a total weight of approx. 5.9 kg (13.01 lbs).

Lift

1

Angled arm enables comfortable carrying

Ergonomic grip for easy carrying

- The arm is angled to keep user's wrists in a comfortable, neutral position.
- The eyepiece can be fixed in place to help prevent damage or loss during transportation.
- A rackless stage and stage cover enable safety and stability during long periods of use.

Carry

2

Fixed eyepieces

Rackless stage and stage cover

Convenient and easily accessible power cable storage

Locking pin for easy binocular rotation

- A locking pin keeps the observation tube in place.
- Interpupillary adjustments, ranging from 48 to 75 mm, enables individual users to set it to their needs.
- Eyepoint adjustments accommodate the user's height and provide greater comfort.
- Left and right diopter adjustment enable optimal focus for each eye.

Prepare

3

Interpupillary adjustment

Eyepoint: height adjustment

## Designed for the Classroom


### Set the Sample

- An inward-facing rotating nosepiece facilitates a larger working area above the stage, enabling easier operation and specimen exchange or making it easier to apply oil to immersion objectives.
- Objectives with long working distances provide clearance for thick samples and slides.
- Bring specimens into focus quickly by using the coaxial knobs for fine or coarse adjustments. The microscope is built for durability and precise control, whether accessed from the left or right side.
- A torque release function helps prevent damage if excess force is applied to the coarse focusing knob or stage knob.

Inward-facing revolving nosepiece


Coaxial knobs

### Observe

- The focus lock helps prevent the objectives from damaging the specimen.
- Field number (FN) of 20 provides a wide field of view.
- The CX23 microscope features plan achromat objectives that preserve outstanding image flatness throughout the field of view.
- The energy saving LED light source provides a long operation lifetime of 20,000 hours.
- Daylight balanced LED illumination preserves the vivid colors on HE-stained samples.

Focus lock

Wide field of view

### Store

- A storage compartment on the back of the CX23 microscope makes it easy to stow the power cable after use.
- Protect your investment with an optional custom designed wooden case.
- When the CX23 microscope needs to be left unattended, a built-in security slot enables attachment of an anti-theft cable.

Cable storage compartment


Built-in security slot

## CX23 Specifications

Model		CX23LEDRF (S <sup>+</sup> /S <sup>2</sup> )	CX23RTFS2
Body		Aluminum die-cast metal frame, protective plastic covering	
Optical System		Infinity optical system	
Illumination System		Built-in transmitted illumination system LED power consumption: 0.5 W (nominal value)	
Focusing		<ul style="list-style-type: none"> <li>• Stage height movement</li> <li>• Coarse movement stroke: 15 mm</li> <li>• Coarse adjustment limit stopper</li> <li>• Torque adjustment for coarse adjustment knob</li> <li>• Fine focus knob (minimum adjustment gradations: 2.5 μm)</li> </ul>	
Revolving Nosepiece		Fixed quadruple nosepiece with inward tilt	
Stage		Wire movement mechanical fixed stage <ul style="list-style-type: none"> <li>• Size (W × D): 174 mm × 89 mm</li> <li>• Travelling range (X × Y): 76 mm × 30 mm</li> <li>• Specimen holder</li> <li>• Specimen position scale</li> </ul>	
Observation Tube	Type	Binocular, anti-fungal	Trinocular, anti-fungal
	Eyepiece (anti-fungal)	10X Field Number (FN): 20 15X Eyepiece (FN): *2 (optional)	
	Tube inclination	30°	
	Light path selector	None	Eyepiece/Camera port = 100/0 or 0/100
	Interpupillary distance adjusting range	48 – 75 mm	
Eyepoint adjustment		370.0 – 432.9 mm	
Condenser		• Abbe condenser NA 1.25 with oil immersion	• Built-in aperture iris diaphragm
Objectives		Plan Achromat, anti-fungal 4X NA 0.1 W.D. 27.8 mm 10X NA 0.25 W.D. 8.0 mm 20X NA 0.4 W.D. 2.5 mm (optional) 40X NA 0.65 W.D. 0.6 mm 100X Oil NA 1.25 W.D. 0.13 mm (CX23LEDRFS1 built-in, CX23LEDRFS2/RTFS2 optional)	
Optional Accessories		Reflection mirror (CH20-MM), eyepiece micrometer Dust cover, darkfield stop (CH2-DS+CH2-FH) Field Diaphragm (CX23-FS)* Dedicated wooden case (CH20-WB2)	Reflection mirror (CH20-MM), eyepiece micrometer Dust cover, darkfield filter (CH2-DS)
Dimensions (W × D × H)		198 mm × 398 mm × 386 mm	198 mm × 398 mm × 430 mm
Weight		Approximately 5.9 kg	Approximately 6.5 kg
Rated Voltage/Electric Current		AC 100 – 240 V 50/60 Hz 0.4 A	
Power Consumption		Less than 2 W	

\*CX23-FS is included in CX23RTFS2 standard set

## DIMENSIONS

Dedicated wooden case (optional)

Specimen courtesy of :  
 JAPANESE FOUNDATION FOR CANCER RESEARCH  
 Cancer Institute, Cancer Institute Hospital  
 Department of Pathology  
 Futoshi Akiyama, M.D., Ph.D.  
 Yuichi Ishikawa, M.D., Ph.D.

- OLYMPUS CORPORATION is ISO14001 certified.
- OLYMPUS CORPORATION is ISO9001 certified.
- OLYMPUS CORPORATION is ISO13485 certified.

\* All company and product names are registered trademarks and/or trademarks of their respective owners.  
 • Specifications and appearances are subject to change without any notice or obligation on the part of the manufacturer.

[www.olympus-lifescience.com](http://www.olympus-lifescience.com)

**OLYMPUS**<sup>®</sup>

OLYMPUS CORPORATION  
 Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-3914, Japan

Printed in Japan N8600083-032017

**Motic**  
MORE MICROSCOPY

## Moticam 1080 HDMI & USB

HD CAMERAS

The HD affordable multi-tasking microscopy camera. The Moticam 1080 is dedicated to documentation of microscopic results with a clear focus on maximum colour fidelity and fast live image.

It does not require a computer for standard operation, as it comes with an on-board software that will allow you to use it by simply connecting it to a screen through its HDMI port, and with a mouse connected to the USB port of the camera.

[www.moticamseries.com](http://www.moticamseries.com) | [www.moticeurope.com](http://www.moticeurope.com)

### FULL HD STANDALONE CAMERA

## Moticam 1080 TECHNICAL SPECIFICATIONS


<b>Sensor Type</b>	CMOS
<b>Sensor Size</b>	1/2.8"
<b>Resolution</b>	8.0 MP
<b>Capture format (on SD-card)</b>	Still image 3840 x 2160(8.0MP) / Video Full HD 1980 x 1080
<b>Live Display Mode (through USB)</b>	1920 x 1080 (Full HD) @ 30 fps*
<b>Live Display Mode (through HDMI)</b>	1920 x 1080 (Full HD) @ 60 fps*
<b>Pixel Size</b>	2.8µm x 2.8µm
<b>Data Transfer</b>	HDMI (1080p) and USB2.0
<b>Slot</b>	SD Card (max 32Gb)
<b>Supported OS</b>	Microsoft Windows XP SP3/Vista/7/8/10 and MAC OSX & Linux
<b>Minimum Computer Requirements</b>	2GHz dualcore – RAM memory 2GB – Video memory min. 512 MB
<b>Lens Mount</b>	C-Mount
<b>On-board software</b>	Standard tools
(Camera connected through HDMI to a HD screen and through USB to a mouse)	Still image capture, Live image record, Zoom in & out, Magnifier, ROI, Mirror, Rotation, Split, Freeze, Cross line, Gallery
	<b>Image adjustments (AE/AWB)</b>
	Auto/Manual Exposure mode, Gain, Target Brightness setting, Light Frequency, Auto/Manual White Balance mode.
	<b>Advanced settings</b>
	Saturation, Contrast, Gamma, Sharpness, Denoise, User settings, System's Time setting, Language (English, Chinese and Japanese), Format SD card.
<b>Package Includes</b>	CS Ring Adaptor, HDMI cable, Focusable Lens, 30mm and 38mm Eyepiece Adapters, USB 2.0 Cable, Calibration Slide, Macro Tube, Motic Images Plus 3.0 for PC/OSX/Linux

\*frames per second under optimal illumination conditions

### ON-BOARD SOFTWARE & MOTIC IMAGES PLUS 3.0 SOFTWARE

#### On-board software

The Moticam 1080 has its own built-in software, making it possible to view live images, make image adjustments or edit them. You can also capture images or record videos directly from a screen with the help of a mouse, storing them in an SD card.

#### Motic Images Plus 3.0

By connecting your Moticam to your computer you will be able to work with our well-known software, that comes included in the package. View, capture, edit, measure, make reports... all its standard features and the new ones packed in a new user-friendly interface.

# Moticam 1080 HDMI & USB

<b>Sensor Type</b>	CMOS	
<b>Sensor Size</b>	1/2.8"	
<b>Resolution</b>	8.0 MP	
<b>Imaging Area</b>	6.4 mm (Diagonal)	
<b>Capture format (on SD-card)</b>	Still Image 3840 x 2160 (8.0MP) / Video Full HD 1920 x 1080	
<b>Live Display Mode (through USB)</b>	1920 x 1080 (Full HD) @ 30 fps*	
<b>Live Display Mode (through HDMI)</b>	1920 x 1080 (Full HD) @ 60 fps*	
<b>Pixel Size</b>	2.8µm x 2.8µm	
<b>Scan Mode</b>	Progressive	
<b>Data Transfer</b>	HDMI (1080p) and USB2.0	
<b>Exposure Time</b>	0.06ms ~ 40ms with Auto Exposure mode 0.06ms ~ 999ms with Manual Exposure mode	
<b>Operating Temperature</b>	From -10 to +60 Degrees Celsius non condensing	
<b>Max. Signal to Noise Ratio</b>	40dB	
<b>Dynamic Range</b>	71dB	
<b>Sensitivity</b>	0.51V (Typical at 1/30 sec)	
<b>Slot</b>	SD Card (max 32Gb)	
<b>Buttons</b>	On / Off	
<b>Support Device</b>	Twain and Direct Show Driver	
<b>Supported OS</b>	Microsoft Windows XP SP3/Vista/7/8/10 and MAC OSX & Linux	
<b>Minimum Computer Requirements</b>	2GHz dualcore – RAM memory 2GB – Video memory min. 512 MB	
<b>Lens Mount</b>	C-Mount	
<b>Focusable Lens</b>	12mm	
<b>On-board software</b>	Standard tools	Still Image capture, Live image review, Zoom in & out, Magnifier, ROI, Mirror, Rotation, Split, Freeze, Cropping, Softkey
(Camera connected through HDMI to a HD screen and through USB to a mouse)	Image adjustments (AE/AWB)	Auto/Manual Exposure mode, Gain, Target Brightness setting, Light Frequency, Auto/Manual White Balance mode
	Advanced settings	Resolution, Contrast, Gamma, Shutterless, Defocus, Base settings, System's Time setting, Language (English, Chinese and Japanese), Format SD Card
<b>Power supply</b>	12V (External Power supply)	
<b>Package Includes</b>	CS Ring Adaptor, HDMI cable, Focusable Lens, 30mm and 38mm Eyepiece Adapters, USB 2.0 Cable, Calibration Slide, Macro Tube, Motic Images Plus 3.0 for PC/OSX/Linux	

\*frames per second under optimal illumination conditions