

**EXCELLENT
ENVIRONMENTAL
FOOTPRINT**

Power consumption of the machine even when exposing is comparable to that of a laptop.

WORLDWIDE USER BASE

Over 170 laboratories around the world, including national labs and international leading Universities.

**INTUITIVE WINDOWS®
USER INTERFACE**

Designed for use by PhD students and post-docs in a research environment while offering high levels of flexibility.


**COMPETITIVE PRICE
AND LOW COST OF
OWNERSHIP**

Affordable price ideal for universities and industrial R&D.

For more information, please visit us on the Web at:

www.durhammagnetooptics.com

MicroWriter ML® 3 family overview



The MicroWriter ML® products are a range of photolithography machines designed for rapid prototyping and small volume manufacturing in R&D laboratories and clean rooms.

Conventional approaches to photolithography are usually based on exposing through a chromium-glass mask manufactured by specialist vendors. In R&D environments it is often necessary to change the mask design frequently. Direct-write lithography tools (also known as digital mask aligners or maskless aligners) overcome this problem by holding the mask in software. Rather than projecting light through a physical mask, direct-write lithography uses computer-controlled optics to project the exposure pattern directly onto the photoresist.

The MicroWriter ML®3 family comprises four compact, high-performance, direct-write optical lithography machines which are designed to offer unprecedented value for money in a small laboratory footprint. All machines have an excellent environmental footprint: power consumption of the machine even when exposing is comparable to that of a laptop.

DURHAM MAGNETO OPTICS LTD
BUILDING 7, OLD FARM BUSINESS CENTRE, CHURCH ROAD, TOFT, CAMBRIDGE CB23 2RF, UK.

DURHAM MAGNETO OPTICS INC
2530 MERIDIAN PARKWAY, 3RD FLOOR, DURHAM, NORTH CAROLINA 27713, USA.



FAST WRITING SPEEDS

120mm²/minute at 2µm resolution and 180mm²/minute at 5µm resolution, allowing a typical 100mm x 100mm area to be exposed at 2µm resolution in under 2 hours.

PROPER AUTOMATIC LENS CHANGER

Automatically changes microscope objective lenses and exposure resolution beams using a motorised motor.

DUAL WAVELENGTH EXPOSURE

LIGHTSOURCE OPTION

Adds both 365nm light source and 405nm light source; software selectable.

BACKSIDE ALIGNMENT OPTION

Adds imaging module underneath the wafer chuck, allowing alignments of structures on both sides of wafers.

For more information, please visit us on the Web at:
www.durhammagnetooptics.com

MicroWriter ML® 3 Baby

This is our entry level machine and is one of the lowest cost direct-write optical lithography machines available anywhere in the world. It operates at a single minimum feature size of 1µm with a wavelength of 405nm and is designed to sit on a standard laboratory bench either in a clean-room or in a general laboratory. A high quality optical microscope with a x10 Olympus objective allows exposures to be aligned to existing structures or to the edges of the substrate. Despite its low cost, it is still fast with a top writing speed of 50mm²/minute, allowing a typical 50mm x 50mm area to be exposed in under 1 hour.

MicroWriter ML® 3 Baby Plus

The Baby Plus adds a number of features to the Baby which are usually only found in high-end machines. Two different minimum feature sizes (1µm and 5µm) can be selected automatically via software and without the user needing to exchange any lens manually. This allows non-critical parts of the exposure to be performed rapidly at 5µm minimum feature size while retaining high resolution writing for critical parts. Locating alignment markers or edges of substrates is faster thanks to an automatic lens changer on the optical microscope allowing the user to switch between x3 and x10 objectives via software. It also features an optical surface profilometer tool and an automated wafer inspection tool for examining fabricated structures. An edge locator allows wafers and dies to be centred automatically. Writing speeds are some of the fastest on the market: up to 50mm²/minute at 1µm minimum feature size and up to 180mm²/minute at 5µm minimum feature size, allowing a typical 50mm x 50mm area combining critical and non-critical areas to be exposed in under 30 minutes.

MicroWriter ML® 3 Mesa

The Mesa has all of the features of the Baby Plus and adds a 0.6µm minimum feature size lens and x20 microscope objective, making it a table-top lithography tool with sub-micron resolution.

MicroWriter ML® 3 Pro

This is our flagship machine and best seller and offers no-compromise sub-micron lithography on up to 9" wafers. It is designed for highly demanding individual research groups or for central clean room facilities. Four different minimum feature sizes (0.6µm, 1µm, 2µm and 5µm) can be selected automatically via software and without the user needing to exchange any lens manually. The optical microscope contains a full set of high performance bright and sharp infinite conjugate objectives (x3, x5, x10 and x20) with a software controlled automatic lens changer, allowing large substrate areas to be searched rapidly and individual sub-micron objects such as nanowires and crystal flakes to be accurately located. An additional lens offering 0.4µm minimum feature size and x50 microscope is available as an option. Top writing speeds are very fast: 15mm²/minute at 0.6µm resolution, 50mm²/minute at 1µm resolution, 120mm²/minute at 2µm resolution and 180mm²/minute at 5µm resolution. In addition to the optical surface profilometer tool and automated wafer inspection tool present in the Baby Plus and Mesa, there is also a Virtual Mask Aligner mode in which the pattern to be exposed is displayed on top of the real-time microscope image, allowing the machine to be used like a traditional mask aligner. A backside alignment camera for aligning double-polished wafers is available as an option.



AUTOCALIBRATION

Allows users to check and correct calibration.

AUTOMATIC MARKER RECOGNITION

Automatically identifies the precise position of lithographic markers visible the real-time microscope.

AUTOMATIC BARCODE GENERATION AND RECOGNITION

Automatically creates the exposure pattern for 2D barcode. Developed barcode can be identified automatically.

COMPACT LABORATORY FOOTPRINT

70cm (w) x 70cm (d) x 75cm (h) for a Baby, a Baby Plus, or a Mesa.

90cm (w) x 75cm (d) x 153cm (h) (including optical table; excluding PC workstation) for a Pro.

TECHNICAL SUPPORT


International network of trained local service engineers to keep you running.

For more information, please visit us on the Web at:


www.durhammagnetooptics.com

Why choose the MicroWriter ML[®] 3 family?


Virtual mask aligner.



Combining resolutions.



Greyscale lithography.



Wide field viewer.



Barcode generation.



Multiple samples.




Wafer inspection.



Producing photomasks.



Structures with a high aspect ratio.



Well-established user base around the world.



- All of our machines are very competitively priced.
- All of our machines have fast writing speeds.
- All of our machines have a low cost of ownership. Our lightsources have a lifetime of 20,000 hours and are guaranteed for 5 years.
- The MicroWriter ML[®]3 Baby Plus, MicroWriter ML[®]3 Mesa and MicroWriter ML[®]3 Pro have an impressive array of advanced features usually only found in high-end machines.
- All of our machines are designed for use by PhD students and post-docs in a research environment and so have an attractive, intuitive and simple Windows[®] user interface while offering the flexibility and high levels of access to machine operation for those who want to develop new techniques.
- All of our machines handle the small millimetre-size chips and delicate objects often used in R&D, as well as large wafers.
- All of our machines provide optical autofocusing as standard. This works well for large wafers as well as for small dies and other unusual substrates.
- All of our machines have a high performance laser interferometer as standard.
- Proper dual wavelength option gives users the best flexibility between high resolution and SU-8 patterning.
- All of our machines share a common technology platform, allowing you to upgrade from MicroWriter ML[®]3 Baby to MicroWriter ML[®]3 Baby Plus, MicroWriter ML[®]3 Mesa and to MicroWriter ML[®]3 Pro at a later date.
- There is a well-established user base of MicroWriter ML[®] machines in over 170 laboratories around the world, including national labs and internationally leading Universities.
- We have an international network of trained local service engineers to keep you running.

Durham Magneto Optics



DURHAM MAGNETO OPTICS LTD
BUILDING 7, OLD FARM BUSINESS CENTRE, CHURCH ROAD, TOFT, CAMBRIDGE CB23 2RF, UK.

DURHAM MAGNETO OPTICS INC
2530 MERIDIAN PARKWAY, 3RD FLOOR, DURHAM, NORTH CAROLINA 27713, USA.

ALL OPTICAL AUTOFOCUS

Works well on large wafers as well as on small samples or on non-conventional samples.

AUTOMATIC WAFER CENTERING

Allows exposures to be centred on the wafer automatically.

EXPORT IMAGE TOOL ("DRAW MODE")

Allows exposures to be designed directly on top of an image taken from the real-time microscope.

WAFER INSPECTION TOOL

Automatically visits a list of user-supplied coordinates and takes a microscope image.

ACCESSING SYSTEM USAGE TIME

Stores usage time data of different users' accounts.

For more information, please visit us on the Web at:
www.durhammagnetooptics.com

Detailed comparison of features and performance

	MicroWriter ML® 3 Baby	MicroWriter ML® 3 Baby Plus	MicroWriter ML® 3 Mesa	MicroWriter ML® 3 Pro
Maximum substrate size.	155mm x 155mm x 7mm.	155mm x 155mm x 7mm.	155mm x 155mm x 7mm.	230mm x 230mm x 15mm. 330mm x 330mm x 175mm as an option.
Maximum writing area.	149mm x 149mm.	149mm x 149mm.	149mm x 149mm.	195mm x 195mm. 295mm x 295mm as an option.
Exposure minimum feature sizes.	1µm.	1µm and 5µm.	0.6µm, 1µm, and 5µm.	0.6µm, 1µm, 2µm, and 5µm. 0.4µm as an option.
Surface tracking autofocus system?	Yes.	Yes.	Yes.	Yes.
Edge locating module for automatic wafer centring?	Yes.	Yes.	Yes.	Yes.
Greyscale lithography?	Yes, up to 255 grey levels.	Yes, up to 255 grey levels	Yes, up to 255 grey levels.	Yes, up to 768 grey levels.
Export image tool ("Draw mode")?	Yes.	Yes.	Yes.	Yes.
Alignment microscope objectives.	x10.	x3 and x10.	x3, x10, and x20.	x3, x5, x10, and x20. x50 as an option.
Automatic lens changer for exposure resolution and alignment microscope?	No.	Yes.	Yes.	Yes.
Backside alignment?	No.	No.	No.	As an option.
Exposure wavelength.	405nm. 385nm, 365nm, and dual wavelength (405nm-365nm) as an option.	405nm. 385nm, 365nm, and dual wavelength (405nm-365nm) as an option.	405nm. 385nm, 365nm, and dual wavelength (405nm-365nm) as an option.	385nm. 405nm, 365nm, and dual wavelength (405nm-365nm) as an option.
Maximum writing speed.	50mm ² /minute at 1µm resolution.	50mm ² /minute at 1µm resolution and 180mm ² /minute at 5µm resolution.	15mm ² /minute at 0.6µm resolution, 50mm ² /minute at 1µm resolution, and 180mm ² /minute at 5µm resolution.	15mm ² /minute at 0.6µm resolution, 50mm ² /minute at 1µm resolution, 120 mm ² /minute at 2µm resolution, and 180mm ² /minute at 5µm resolution.



OPTICAL SURFACE PROFILOMETER TOOL

Allows the user to measure small variations in height of a sample.

VIRTUAL MASK ALIGNMENT OPTION

Allows the user to see precisely where the pattern will be exposed.

MULTIPLE WAFERS/CHIPS OPTION

Different users can run their own exposures automatically overnight.

GLOVE BOX OPTION

Run exposures under an inert gas environment.

COMMON TECHNOLOGY PLATFORM

Allows the users to upgrade the MicroWriter at a later date.

For more information, please visit us on the Web at:
www.durhammagnetooptics.com

Detailed comparison of features and performance

	MicroWriter ML® 3 Baby	MicroWriter ML® 3 Baby Plus	MicroWriter ML® 3 Mesa	MicroWriter ML® 3 Pro
Overlay alignment accuracy at best resolution.	±2µm.	±1µm.	±0.5µm.	±0.5µm.
Minimum addressable grid.	100nm.	60nm.	30nm.	30nm. 12.5nm as an option.
Motion stage minimum XY step size.	15nm.	15nm.	15nm.	4nm.
XY interferometer resolution.	15nm.	15nm.	15nm.	1nm.
Optical surface profiler Z resolution.	Not applicable.	200nm.	200nm.	100nm.
Autocalibration tool?	Yes.	Yes.	Yes.	Yes.
Automatic marker recognition tool?	Yes.	Yes.	Yes.	Yes.
Automatic wafer inspection tool?	No.	Yes.	Yes.	Yes.
Virtual Mask Aligner tool?	No.	As an option.	As an option.	Yes.
Temperature compensated enclosure?	No.	As an option.	As an option.	Yes.
Can handle multiple wafer/chips?	No.	As an option.	As an option.	Yes.
Can be installed inside a glove box?	As an option.	As an option.	As an option.	As an option.
Supplied with vibration isolating optical table?	No.	No.	No.	Yes.
Mask design software?	Open source KLayout supplied. Clewin as an option.	Open source KLayout supplied. Clewin as an option.	Open source KLayout supplied. Clewin as an option.	Clewin supplied.
Can be upgraded to higher members of the MicroWriter ML® 3 family?	Yes.	Yes.	Yes.	Not applicable.



FRICTION CHUCK

Carefully designed friction chuck allows MEMS devices with nitride windows or other delicate substrates to be used; no minimum wafer size.

CURVED SUBSTRATES

Perform exposures across a variety of substrates, including flat and curved forms, Si, glass, ceramic, diamond, and liquid polymers.

PHOTOMASKS

Produce photomasks conveniently and cheaply.

FREE SOFTWARE UPGRADE

Receive free software upgrades for the lifetime of the machine.


COMPANY CULTURE AND PHILOSOPHY

We are from a research and design (R&D) background based in Cambridge, UK and the Research Triangle Park, Durham, NC, USA.

For more information, please visit us on the Web at:

www.durhammagnetooptics.com

Examples of fabricated structures



- Scanning Electron Microscope images of resolution-limited structures.
Left: Dots array with diameter of $0.4\mu\text{m}$.
Right: Lines array with width of $0.6\mu\text{m}$.
Structures were produced on Si/LOR/S1805 ($0.5\mu\text{m}$).
- Scanning Electron Microscope images of micro-moulds.
Structures were produced on a $50\mu\text{m}$ thick SU8 layer.
Aspect ratio of the dots array (bottom left) is 8.
- Optical Microscope images of patterns produced across varied types of substrates:
Top left: AlN ceramic.
Top right: Liquid polymer.
Bottom left: Si/SiO₂.
Bottom middle: Glass.
Bottom right: Diamond.
- Optical Microscope images of patterns produced on a photomask.

DURHAM MAGNETO OPTICS LTD
BUILDING 7, OLD FARM BUSINESS CENTRE, CHURCH ROAD, TOFT, CAMBRIDGE CB23 2RF, UK.

DURHAM MAGNETO OPTICS INC
2530 MERIDIAN PARKWAY, 3RD FLOOR, DURHAM, NORTH CAROLINA 27713, USA.



POLOS® SPIN150x SPIN COATER

The new POLOS® SPIN150x offers precise, repeatable process control. It is made from natural polypropylene (NPP) or optionally available in chemically resistant PTFE. This new version of the well-proven model comes in a standard and an advanced version. It is ideal for processing a wide range of substrates with a diameter of up to 150 mm or square samples with dimensions 100 x 100 mm.



UNIQUE DESIGN

The unique outer shell and drain design allows easy switching between table-top and in-deck model. The modular design enables the user to upgrade the unit using a wide range of different accessories to ease dispensing and overall handling.

MOTOR HOMING POSITION

The new SPINx-series has the ability to define a motor homing position, allowing for easy integration in robot controlled/automated environments.

EASY CHAMBER ACCESS

The vacuum chuck sits above the edge of the bowl allowing easy access to the wafer/substrate with an end-effector, tweezers or vacuum wand. This is a unique requirement to enable robot handling.

SPECIFICATIONS HARDWARE:

- Liquid filter trap
- Automatic lid, also controllable via foot pedal (advanced version)
- Programmable motor homing position
- Center injection holder for syringe or dispense nozzle
- Lid lock and vacuum sensor for user safety
- Large (detachable) touchscreen display
- USB-port to store recipes on USB-drive and for software updates specifications drive-unit
- Indirect brushless drive unit - up to 12.000 RPM
- High acceleration and accuracy: 1 - 30.000 RPM
- Clockwise/counter clockwise rotation and puddle mode
- Unique design to switch between desktop and in-deck model

AUTO LID

The lid can be automatically opened and closed using the User Interface or alternative a foot pedal (ideal for glovebox usage). In addition, it is also possible to automatically open the lid as a final step of your recipe.

LIQUID FILTER TRAP

The SPINx-series is equipped with a liquid filter trap to protect the spin coater's critical components. It will capture any liquids or resists entering the vacuum lines via the process chamber or vacuum chuck in a liquid container. The container can be viewed through a cut-out in the spinner housing, for easy maintenance.

CONFIGURATION

Process chamber Material	Natural polypropylene (NPP) or PTFE
Max. substrate diameter	Up to 6" (150 mm) wafers Up to 4" x 4" (100 mm) substrates

SUITABLE FOR:

- Coating
- Cleaning
- Rinse/Dry
- Developing
- Etching
- PDMI and other processes

STANDARD ACCESSORIES:



- Vacuum chuck for 4 - 6-inch wafers
- Fragment adapter for fragments of 10 mm and larger

SPIN COATING

STANDARD VS. ADVANCED VERSION

MODEL	SPIN150x STANDARD	SPIN150x ADVANCED
Liquid filter trap	Yes	Yes
Unique outer shell and drain design	Yes	Yes
Programmable motor home position	Yes	Yes
Center injection holder for syringe or dispense nozzle	Yes	Yes
Lid lock and vacuum sensor for user safety	Yes	Yes
Large (detachable) touchscreen display	Yes	Yes
USB port to store recipes on USB drive	Yes	Yes
Dispense vessel	No	Yes
Auto-closing lid	No	Yes (optional)
Linear dispense arm	No	Yes (optional)
Backside rinse	No	Yes (optional)
Available in chemical resistant PTFE	No	Yes (optional)

DIMENSIONS SPIN150x



Unique outer/inner shell and drain design allows switch between table-top and in-deck model



DESKTOP MODEL



IN-DECK MODEL

OPTIONS

Syringe holder
starter kitTransparent 0,5 mm
PET liners

Centering tool

Center dispense
system: opusCorrugated
DrainhosePOLOS® Vacuum
pumpPOLOS® Peristaltic
pumpCentral dispensing
syringe holder

And more. Contact us for all options!

POLOS® HOTPLATE 200

Our all NEW table-top hotplate is a versatile and affordable tool for R&D and pilot lines. It is designed with a soft-close lid and is suitable for soft bake as well as hard bake processes, and curing of photoresist, epoxy or any other work requiring precise temperature control. The POLOS® Hotplate 200 is also available as an advanced version. This unit is standard equipped with complimentary features such as an N2 connector, lifting pins, vacuum bake, and proximity pins.

NEW!




Also available with a
350 x 350 mm heating area!

HOTPLATE 200 STANDARD


- Heater surface area 220 x 220 mm
- Suitable for 1 x 200 mm wafer
- Soft-close lid, enabling single-handed closing
- Glove-friendly operation
- Temperature range 50 - 230°C
- Programmable storage of 20 programs (Temperature/Time)
- Countdown timer (1 - 999 sec.)
- Temperature uniformity $\pm 1^\circ\text{C}$
- Voltage: 230 or 110 VAC
- Heater block material: Aluminum (anodized)
- Housing material: Powder coated steel
- Weight: 12 kg
- Dimensions: approx. 422 x 295 x 201 mm

EASILY CONVERTIBLE INTO AN INDECK MODEL

Using our specially designed indeck bracket and display bracket, you can easily convert the POLOS® hotplate into an indeck solution.




INDECK DISPLAY
BRACKET



INDECK BRACKET

HOTPLATE 200 ADVANCED

- Heater surface area 220 x 220 mm
- Suitable for 1 x 200 mm wafer
- Temperature range 50 - 230°C
- Programmable storage of 20 programs (Temperature/Time)
- Countdown timer (1 - 999 sec.)
- Temperature uniformity $\pm 1^\circ\text{C}$
- Equipped with programmable (electric) lifting pins set in radius of 82 mm
- Equipped with proximity pins to hold the wafer above the heating plate while baking
- Equipped with perforated vacuum plate to realize a hard contact bake
- Voltage: 230 or 110 VAC
- Heater block material: Aluminum (anodized)
- Housing material: Powder coated steel
- Including soft-close lid with N2 connection
- Weight: 12 kg
- Dimensions: approx. 422 x 295 x 201 mm



Perforated vacuum plate and programmable lifting pins

SYSTEM DRAWINGS

