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
## PDControl

Software for particle size and concentration analysis with the welas® digital opt. aerosol spectrometer

## Description

The **PDControl software** is designed for measurement value acquisition and data analysis by all of the welas® digital systems. It enables quick and easy particle size and concentration analysis.

The PDControl software provides special advantages for the analysis of rapid changes in concentration and particle size distributions using the high time resolution of 10 ms, depending on the measurement system. In addition, PDControl enables the transparent display of measured distributions with analysis of the coincidence percentage as per Dr. Umhauer and Prof. Sachweh (**see last figure**). Palas® software solutions in the fields of particle measurement technology, filter testing, and test rig control are based on our many years of experience and have been continuously optimized and further developed in close collaboration with our international customers based on their requirements. As a result, Palas® software is especially practical and provides many advantages in daily applications.

We will be happy to send you our demo software.


**Palas® also offers other software, including:**

- **FTControl** for filter testing
- **In addition, Palas® offers custom software solutions for the automation and control of customer-specific test rigs with the inclusion of external measurement systems.**


## Quality to the last detail

### PDControl software

- Display of cumulative and distributive size distributions in diagram and table form

- Comparison of cumulative and distributive size distributions in diagram form


- Display of all measured single particles and the corresponding statistics for distributions in table form

dN/N

differential distribution | integral distribution | **table** | statistics

Xu [µm]	X [µm]	Xo [µm]	dx [µm]	dN [P]	dN/N	dN/N/dx [1/µm]	dN/N/log(dx)	Sum(dN) [P]
0,178	0,184	0,191	0,013	1,0	0,00018	0,0134	0,0057	1,0
0,191	0,198	0,205	0,014	2,0	0,00035	0,0249	0,0113	3,0
0,205	0,213	0,221	0,015	6,8	0,00120	0,0782	0,0383	9,8
0,221	0,229	0,237	0,016	11,4	0,00202	0,1227	0,0646	21,1
0,237	0,246	0,255	0,018	11,0	0,00196	0,1108	0,0627	32,2
0,255	0,264	0,274	0,019	16,6	0,00294	0,1546	0,0940	48,8
0,274	0,284	0,294	0,020	16,0	0,00284	0,1392	0,0910	64,8
0,294	0,305	0,316	0,022	13,3	0,00235	0,1071	0,0753	78,1
0,316	0,328	0,340	0,024	14,2	0,00252	0,1060	0,0807	92,3

- Storage/loading of measurements
- Read-only mode for measurement files
- Data export in a text file (possible to import data into Excel)
- Acquisition and analysis in as little as 10 ms ideal for rapid, processes
- Up to four statistical values are able to be displayed simultaneously, e.g. number concentration, coincidence percentage, Sauter diameter, and mean diameter


### Coincidence and border zone error analysis

welas® digital is equipped with a new, fast 20 MHz signal processing processor that analyzes the curve of each individual particle signal. This enables coincidences, i.e. more than one particle at the same time in the measurement volume, to be detected on the individual signal and corrected (as per Dr. Umhauer / Prof. Sachweh).

The PDControl software enables our customers to trace each measurement detail, using:

- Display of measured distributions in diagram form with
  - Valid particle size distribution
  - Border zone and coincidence distributions and
  - Statistically corrected particle size distributions.

## Benefits

- Optimal information content:
  - Display of distributions in diagram and table form: Number, area, volume, concentration (number and mass)
  - Comparison of multiple distributions in one diagram
  - Display of the time sequence of statistical values
  - Analysis of 14 statistical values
- Acquisition and analysis in as little as 10 ms ideal for rapid, unsteady processes
- Coincidence analysis
- Simultaneous measurement with multiple welas® digital systems
- Also able to be used with other particle measuring devices
- Acquisition by four external sensor signals in parallel with particle measurement, e.g. differential pressure, temperature, relative humidity, volume flow, etc.
- Clear structure
- Easy installation on the PC
- Short processing times
- Regular updates via the Internet free of charge
- Easy to transfer data to Excel
- Unlimited copies for various workstations within a company
- Easy to operate
- Reliable function
- Reduces your operating expenses

## Datasheet

<i>Parameter</i>	<i>Description</i>
<b>Operating system</b>	Windows® 2000/ Windows® XP, Windows® Vista, Windows® 7
<b>Processor</b>	Min. Pentium 4, 1800 MHz
<b>User memory</b>	From 4 GB
<b>Screen resolution</b>	At least 1280 • 800

## Applications

- Particle measurement and particle size analysis
- Measurement of MDI/DPI inhalers and nebulizers
- Emission and immission measurements
- Analysis of rapid, unsteady processes
- Smoke detector tests
- Particle measurement for clouding processes

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
## Promo<sup>®</sup> 2000


Scattered-light aerosol spectrometer system for process measurement technology and monitoring applications

## Model Variations


### Promo<sup>®</sup> 2000 H

With heating regulation up to 250 °C for welas<sup>®</sup> aerosol sensors

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## Description

Promo® 2000 is a light-scattering aerosol spectrometer system for particle size analysis and concentration determination that can be equipped with all *welas*® sensors.

On Promo® 2000, the *welas*® sensors equipped with different measurement volumes, as required, can be easily connected via fiber optic cables and interchanged as required. These sensors allow reliable measurement in the concentration range from  $< 1$  particle/cm<sup>3</sup> up to  $10^6$  particles/cm<sup>3</sup> and are available for measurement in gases, as well as in liquids. Please observe the "*welas*® sensors" product data sheet.

Unique are up to four measuring ranges in only one device:

- 0.2 µm – 10 µm
- 0.3 µm – 17 µm
- 0.6 µm – 40 µm
- 2 µm – 100 µm (additionally for sensors 2300 and 2500).

Promo® 2000 is famous for up to 128 size channels per measuring range and a concentration range from  $< 1$  particle/cm<sup>3</sup> to  $10^6$  particles/cm<sup>3</sup>.

A touch display ensures user-friendly operation. Measurements can be started easily, and all data, such as the current number distribution and the number concentration, as well as 24 further statistical values, can be evaluated and displayed in real time.

With Promo® as a standalone measuring device (i.e. without external computer) measurements are performed continuously. All incoming data can be stored with a maximum temporal resolution of 1 s. Promo® 2000 can therefore measure and save data over weeks independently. For data transfer, Promo® can also be integrated into a company network.

Promo® 2000 has a standard interface and can be controlled by a process control system or by a simple Labview program. For this reason, Promo® 2000 **is especially well suited for control and monitoring applications**. Temperature, humidity and pressure sensors can be connected.

Palas® offers remote maintenance for the device and data access via [www.palas.de/user](http://www.palas.de/user).


### The new Promo® 2000 measurement technology:

Promo® 2000 offers a new, fast 20 MHz signal processing processor, which analyses the progression of each particle signal. This makes it possible to recognise coincidental events in light scattering measurement technology at the individual signal and correct them (according to Dr. Umhauer / Prof. Dr. Sachweh). This way it is possible to increase the maximum concentration limit up to  $10^6$  particles/cm<sup>3</sup> (*welas*®

2070 sensor). Also in low concentrations  $< 1$  particle/cm<sup>3</sup> with the welas<sup>®</sup> 2500 sensor, this leads to a higher measuring accuracy.


**High classification accuracy and high particle size resolution are guaranteed by the following special features (see Graph 1):**

- White light and 90° light-scattering detection  
⇒ **Unambiguous calibration curve**
- Patented T-aperture  
⇒ **No border zone error**
- New digital individual signal processing  
⇒ **Coincidence detection and correction of the individual signal making it possible to measure higher concentrations**

Graph 1: Example with 2200 sensor

Promo<sup>®</sup> 2000 is characterized by its very high counting efficiency starting from 0.2 μm!

Graph 2: Example with 2200 sensor, in relation to LAS-X II

**The new Promo<sup>®</sup> 2000:**

The Promo<sup>®</sup> 2000 is characterized by its optical fibre technology. The welas<sup>®</sup> aerosol sensor is connected via a fiber optic cable with a length up to 30 m with the Promo control unit. This leads to a minimization of particle losses in long sampling lines by simply installing the sensor directly at the sampling location.

## Benefits

- Measuring range of 0.2 to 100 µm (4 measuring ranges selectable in one device)
- Up to four measuring ranges in only one device:
  - 0,2 µm – 10 µm
  - 0,3 µm – 17 µm
  - 0,6 µm – 40 µm
  - 2 µm – 100 µm (additionally for sensors 2300 and 2500)
- Up to 128 size channels per measuring range
- Concentration range of 1 particle/cm<sup>3</sup> to 10<sup>6</sup> particles/cm<sup>3</sup>
- Calibration curves for different refractive indices
- Very high and reproducible counting efficiency rate starting at 0.2 µm
- **Pressure-resistant up to 10 barg (optional)**
- **Heatable to 250 °C (optional)**
- Optical fibre technology
- Simple operation with a large touch display
- Calibration, cleaning and lamp replacement can all be performed independently by the customer
- External control by RS 232 or Ethernet
- With analysis software PDAnalyze
- Optional: Software PDControl for operation as welas<sup>®</sup> digital available
  
- Low maintenance
- Reliable function
- Reduces your operating expenses

## Datasheet

<i>Parameter</i>	<i>Description</i>
<b>Interfaces</b>	USB, Ethernet, RS232/485, Wi-Fi
<b>Measurement range (size)</b>	0,2 µm – 10 µm 0,3 µm – 17 µm 0,6 µm – 40 µm 2 µm – 100 µm
<b>Size channels</b>	Up to 128 (64/decade)
<b>Measuring principle</b>	Optical light-scattering
<b>Measurement range (number C<sub>N</sub>)</b>	< 1 • 10 <sup>6</sup> Partikel/cm <sup>3</sup>
<b>Time resolution</b>	Up to 1 s
<b>Thermodynamic conditions</b>	10 – 40 °C, -100 – 50 mbarg
<b>Volume flow</b>	5 l/min
<b>Data acquisition</b>	20 MHz processor, 256 raw data channels, digital
<b>Light source</b>	Xenon Bogenlampe 35 W 100 W
<b>Power consumption</b>	35 W 100 W
<b>User interface</b>	Touchscreen, 800 • 480 Pixel, 7"
<b>Power supply</b>	115 – 230 V, 50 – 60 Hz
<b>Housing</b>	Table housing, optionally with mounting brackets for rack-mounting
<b>Dimensions</b>	185 • 450 • 315 mm (H • W • D) (19")
<b>Support options</b>	direct remote access, Palas <sup>®</sup> webserver service
<b>Weight</b>	Approx. 8 kg (control unit), approx. 2.8 kg (sensor)
<b>Operating system</b>	Windows embedded
<b>Data logger storage</b>	4 GB Compact Flash
<b>Software</b>	PDControl, FTControl, PDAnalyze
<b>Installation conditions</b>	+5 – +40 °C (control unit)

## Applications

- Emission monitoring of installations
- Control of grinding and classification processes
- Monitoring of production processes in the food, pharmaceuticals and chemicals industries
- Testing of complete filters, inertial and wet separators or electrostatic precipitators

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