
Idylla™

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- Dietel M. et al. A multicenter validation study of the Idylla™ BRAF Mutation Test on FFPE tissue of malignant melanoma. Poster AACR 2014, Poster AMP 2014.
- Devogelaere B. et al. BRAF V600 mutation testing on FFPE samples using a novel fully integrated molecular diagnostics platform. Poster AACR 2013.
- Janku F. et al. BRAF Mutation Testing with a Novel, Rapid, Fully-Automated Molecular Diagnostics Prototype Platform. Poster AACR 2013.
- Claes B. et al. A novel fully integrated molecular diagnostics platform with broad sample prep capabilities. Poster Sample Prep 2013.
- Vandenbroucke I. et al. Innovative, rapid and easy to use diagnostic multiplex platforms that enable individualized care. Poster CONFO 2013.
- Vandenbroucke I. et al. A rapid and fully automated multiplex assay for KRAS-BRAF mutations with high mutation sensitivity using novel selective amplification and detection technologies. Poster AACR 2014.
- Dietel M. et al. A multicenter validation study of the Idylla™ BRAF Mutation Test on FFPE tissue of malignant melanoma. Poster AACR 2014, Poster AMP 2014.
- Van der Auwera M. et al. Evaluation of a rapid, sensitive and fully automated Idylla™ BRAF Mutation Test starting directly from FFPE samples. Poster BWP 2014.
- Micalessi I. et al. Evaluation of the Idylla™ BRAF Mutation Test performance using the novel fully-automated Idylla™ System. Poster BWP 2014.
- Janku F. et al. BRAF mutation testing with a rapid, fully integrated molecular. Oncotarget 2015.
- Melchior L. et al. Multi-center evaluation of the novel fully-automated PCR-based Idylla™ BRAF Mutation Test on FFPE tissue of malignant melanoma. J Exp Mol Path 2015.
- Micalessi I. et al. Evaluation of the Idylla™ BRAF Mutation Test performance using the novel fully-automated Idylla™ System on FFPE melanoma samples. Poster EMMO 2015.
- Schiefer A-I. et al. Multicenter Evaluation of a Novel Automated Rapid Detection System of BRAF Status in Formalin-Fixed, Paraffin-Embedded Tissues. J of Molecular Diagnostics 2016.
- Alexandre H. et al. Detection of BRAF Mutations Using a Fully Automated Platform and Comparison with High Resolution Melting, Real-Time Allele Specific Amplification, Immunohistochemistry and Next Generation Sequencing Assays, for Patients with Metastatic Melanoma. PLOS ONE 2016.
- Riveiro-Falkenbach E. et al. BRAF V600 mutation detection in tricky melanoma samples using the new Idylla™ BRAF Mutation Test. Poster ESP 2016.

Research applications

- Schiefer A-I. et al. Evaluation of a novel fully automated PCR-based technique for detection of BRAF-status in FFPE tissues. Poster AMP 2014.
- Colling R. et al. Automated PCR detection of BRAF mutations in colorectal adenocarcinoma. J Clin Path 2015.
- Micalessi I. et al. Evaluation of the Idylla™ BRAF Mutation Test performance using the novel fully-automated Idylla™ System on FFPE colorectal cancer samples. Poster EMMO 2015.
- Yeo MK. et al. The usefulness of a novel fully automated PCR-based Idylla test for detection of the BRAF V600E mutation in thyroid tissue: comparison with PNA-clamping PCR, real-time PCR and pyrosequencing. J Clin Pathol 2016.


Disclaimer

Biocartis trademark and logo are trademarks belonging to Biocartis and are used and registered in Europe. Idylla is a registered trademark in the United States (US) and other countries. Idylla trademark and logo are used trademarks belonging to Biocartis. Idylla™ platform and Idylla™ BRAF Mutation Test are CE-marked IVDs in Europe. Idylla™ BRAF Mutation Assay is available for Research Use Only, not for use in diagnostic procedures. Idylla™ is not for sale in USA, Canada and some other countries. Please check availability with the local Biocartis sales representative.

Technical sheet Idylla™ BRAF Mutation Test

BRAF CE IVD

The Idylla™ BRAF Mutation Test, performed on the Biocartis Idylla™ System, is an *in vitro* diagnostic Test for the qualitative detection of V600E/E2/D and V600K/R/M mutations in codon 600 of the *BRAF* gene. The Idylla™ BRAF Mutation Test, from sample-to-result, starts with formalin-fixed, paraffin-embedded (FFPE) human tissue from metastatic melanoma to liberate DNA for subsequent real-time PCR amplification and detection.

Features

BRAF mutation detection

Codon 600	BRAF V600E	(c.1799T>A)
	BRAF V600E2	(c.1799_1800delinsAA)
	BRAF V600D	(c.1799_1800delinsAT, c.1799_1800delinsAC)
	BRAF V600K	(c.1798_1799delinsAA)
	BRAF V600R	(c.1798_1799delinsAG)
	BRAF V600M	(c.1798G>A)
	BRAF Wild Type	(c.1799T)

RNaseP (acting as Sample Processing Control)

Specimen requirements

Sample Type	FFPE tissue sections (5 to 10µm)
Neoplastic cells	≥50%, if less macrodissection is required
Tissue area	50-600mm² (5µm) 25-300mm² (10µm)

Performance

Analytical Sensitivity	1% mutant in wild type background
Between Laboratory Reproducibility (240 results at 3 sites)	100% agreement for 3.5% BRAF V600E 100% agreement for 5% BRAF V600K
Between Lot Reproducibility (120 results on 3 lots)	98.3% agreement for 3.5% BRAF V600E 100% agreement for 5% BRAF V600K
Total turnaround time	
Time	90 minutes

98.9% overall concordance		Routine reference methods, including further analysis by ddPCR, Idylla™ retest, NGS						
		Codon 12	Codon 13	Codon 59	Codon 61	Codon 117	Codon 146	No mutations
Idylla™ KRAS Mutation Assay	Codon 12	142						
	Codon 13		38					
	Codon 59			5				
	Codon 61	1			15			1
	Codon 117					4		
	Codon 146						20	
	No mutation	1	2					133
TOTAL		144	40	5	15	4	20	134
								362

Idylla™ KRAS Posters & Publications

- Maertens G. et al. A solution for same-day extended RAS testing. Poster ESMO 2015
- Vandenbroucke I. et al. A rapid and fully automated multiplex assay for KRAS-BRAF mutations with high mutation sensitivity using novel selective amplification and detection technologies. Poster AACR 2014
- Solassor J. et al. Multi-Center Evaluation of the Fully Automated PCR-Based Idylla™ KRAS Mutation Assay for Rapid KRAS Mutation Status Determination on Formalin-Fixed Paraffin-Embedded Tissue of Human Colorectal Cancer. PLOS ONE 2016

Research applications

- Dario de Biase. et al. Fully automated PCR detection of KRAS mutations on pancreatic endoscopic ultrasound fine-needle aspirates. J Clin Path 2016.


Disclaimer

Biocartis trademark and logo are trademarks belonging to Biocartis and are used and registered in Europe. Idylla is a registered trademark in the United States (US) and other countries. Idylla trademark and logo are used trademarks belonging to Biocartis. Idylla™ platform and Idylla™ KRAS Mutation Test are CE-marked IVDs in Europe. Idylla™ KRAS Mutation Assay is available for Research Use Only, not for use in diagnostic procedures. Idylla™ is not for sale in USA, Canada and some other countries. Please check availability with the local Biocartis sales representative.

Technical sheet Idylla™ KRAS Mutation Test

KRAS CE IVD

The Idylla™ KRAS Mutation Test, performed on the Biocartis Idylla™ system, is an *in vitro* diagnostic Test for the qualitative detection of 21 mutations in codons 12, 13, 59, 61, 117 and 146 of the KRAS gene. The Idylla™ KRAS Mutation Test, from sample-to-result, starts with formalin-fixed, paraffin-embedded (FFPE) human tissue from metastatic colorectal cancers to liberate DNA for subsequent real-time PCR amplification and detection.

Features

KRAS mutation detection

Codon 12 (exon 2)	G12C	(c.34C>T)
	G12R	(c.34G>C)
	G12S	(c.34G>A)
	G12A	(c.35G>C)
	G12D	(c.35G>A)
	G12V	(c.35G>T)
Codon 13 (exon 2)	G13D	(c.38G>A)
Codon 59 (exon 3)	A59E	(c.176C>A)
	A59G	(c.176C>G)
	A59T	(c.175C>A)
Codon 61 (exon 3)	Q61K	(c.181C>A, c.180_181delinsAA)
	Q61L	(c.182A>T)
	Q61R	(c.182A>G)
	Q61H	(c.183A>C, c.183A>T)
Codon 117 (exon 4)	K117N	(c.351A>C, c.351A>T)
Codon 146 (exon 4)	A146P	(c.436G>C)
	A146T	(c.436G>A)
	A146V	(c.437C>T)

KRAS Total (acting as Sample Processing Control)

Specimen requirements

Sample Type	FFPE tissue sections (5 to 10µm)
Neoplastic cells	≥10%, if less macrodissection is required
Tissue area	50-600mm ² (5µm)
	25-300mm ² (10µm)


Performance

Analytical Sensitivity	LOD ≤5% for all KRAS mutations
------------------------	--------------------------------

Technical sheet Idylla™ NRAS Mutation Test

NRAS CE IVD

The **Idylla™ NRAS Mutation Test**, performed on the Biocartis Idylla™ system, is an *in vitro* diagnostic Test for the qualitative detection of 18 mutations in codons 12, 13, 59, 61, 117, 146 of the *NRAS* gene. The Idylla™ NRAS Mutation Test, from **sample-to-result**, starts with formalin-fixed paraffin-embedded (FFPE) human tissue from metastatic colorectal cancers to liberate DNA for subsequent real-time PCR amplification and detection.

Features

NRAS mutation detection

Codon 12 (exon 2)	G12C (c.34G>T)
	G12S (c.34G>A)
	G12D (c.35G>A)
	G12A (c.35G>C)
	G12V (c.35G>T)
Codon 13 (exon 2)	G13D (c.38G>A)
	G13V (c.38G>T)
	G13R (c.37G>C)
Codon 59 (exon 3)	A59T (c.175G>A)
	Q61K (c.181C>A)
Codon 61 (exon 3)	Q61L (c.182A>T)
	Q61R (c.182A>G)
	Q61H (c.183A>C; c.183A>T)
Codon 117 (exon 4)	K117N (c.351G>C; c.351G>T)
Codon 146 (exon 4)	A146T (c.436G>A)
	A146V (c.437C>T)
NRAS Total (acting as Sample Processing Control)	
BRAF Total (acting as Sample Processing Control)	

Specimen requirements	
Sample Type	FFPE tissue sections (5 to 10µm)
Neoplastic cells	≥10%, if less macrodissection is required
Tissue area	50-600mm ² (5µm)
	25-300mm ² (10µm)
Performance	
Analytical Sensitivity	LOD ≤5% for most prevalent NRAS mutations
	100% agreement for 10% NRAS G12D
Between Laboratory Reproducibility (480 results at 3 sites)	100% agreement for 10% NRAS G12V
	100% agreement for 10% NRAS Q61K
	100% agreement for 10% NRAS Q61R
Between Lot Reproducibility (240 results on 3 lots)	100% agreement for 10% NRAS G12D
	100% agreement for 10% NRAS G12V
	100% agreement for 10% NRAS Q61K
	100% agreement for 10% NRAS Q61R
Total turnaround time	
Time	120 minutes

Accuracy - Clinical Performance Evaluation

99.6% overall diagnostic agreement for the *NRAS* gene, was obtained during the clinical performance evaluation comparing Idylla™ with a sequencing-based reference method.

99.6% overall concordance

Reference Method


NRAS		G12A	G12V	G12C	G12D	G12S	G13D	Q61H	Q61K	Q61L	Q61R	No mutation	Total
Idylla™ NRAS Mutation Test	G12A/V	1	4										5
	G12C			2									2
	G12D				8							1*	9
	G12S					0							0
	G13D						1						1
	Q61H							1					1
	Q61K								7				7
	Q61L									3			3
	Q61R										6		6
	No mutation											201	201
Total		1	4	2	8	0	1	1	7	3	6	202	235

Technical sheet Idylla™ NRAS-BRAF Mutation Test

NRAS **BRAF** **CE** **IVD**

The Idylla™ NRAS-BRAF Mutation Test, performed on the Biocartis Idylla™ system, is an *in vitro* diagnostic Test for the qualitative detection of 18 mutations in **codons 12, 13, 59, 61, 117, 146** of the *NRAS* gene and 5 mutations in **codon 600** of the *BRAF* gene.

The Idylla™ NRAS-BRAF Mutation Test, from **sample-to-result**, starts with formalin-fixed, paraffin-embedded (FFPE) human tissue from metastatic colorectal cancers to liberate DNA for subsequent real-time PCR amplification and detection.


Features

NRAS mutation detection	
Codon 12 (exon 2)	G12C (c.34G>T)
	G12S (c.34G>A)
	G12D (c.35G>A)
	G12A (c.35G>C)
	G12V (c.35G>T)
Codon 13 (exon 2)	G13D (c.38G>A)
	G13V (c.38G>T)
	G13R (c.37G>C)
Codon 59 (exon 3)	A59T (c.175G>A)
Codon 61 (exon 3)	Q61K (c.181C>A)
	Q61L (c.182A>T)
	Q61R (c.182A>G)
	Q61H (c.183A>C, c.183A>T)
Codon 117 (exon 4)	K117N (c.351G>C, c.351G>T)
Codon 146 (exon 4)	A146T (c.436G>A)
	A146V (c.437C>T)
NRAS Total (acting as Sample Processing Control)	
BRAF mutation detection	
Codon 600	BRAF V600E (c.1799T>A; c.1799_1800delinsAA)
	BRAF V600D (c.1799_1800delinsAC)
	BRAF V600K (c.1798_1799delinsAA)
	BRAF V600R (c.1798_1799delinsAG)
BRAF Total (acting as Sample Processing Control)	

Technical sheet Idylla™ EGFR Mutation Test

EGFR **CE** **IVD**

The Idylla™ EGFR Mutation Test, performed on the Biocartis Idylla™ System, is an *in vitro* diagnostic test for the qualitative detection of **exon 18 (G719A/C/S), exon 21 (L858R, L861Q), exon 20 (T790M, S768I) mutations, exon 19 deletions and exon 20 insertions** in the *EGFR* oncogene. The Idylla™ EGFR Mutation Test, from **sample-to-result**, starts with formalin-fixed, paraffin-embedded (FFPE) human tissue from non-small cell lung cancer (NSCLC) to liberate DNA for subsequent real-time PCR amplification and detection.

Features

EGFR mutation detection

Exon 18	G719A	c.2156G>C
	G719C	c.2155G>T; c.2154_2155delinsTT
	G719S	c.2155G>A
Exon 19	Del9	c.2238_2248delinsCC
		c.2239_2248delinsC
		c.2240_2248del
		c.2239_2247del
		c.2239_2251delinsC
Exon 20	Del12	c.2240_2251del
		c.2235_2249del
		c.2236_2250del
		c.2239_2253del
		c.2240_2254del
		c.2238_2252del
		c.2237_2251del
		c.2235_2252delinsAAT
		c.2237_2252delinsT
		c.2234_2248del
Exon 21	Del15	c.2236_2253delinsCTA
		c.2237_2253delinsTA
		c.2237_2253delinsTC
		c.2235_2251delinsAG
		c.2236_2253delinsCAA
Exon 22	Ins20	c.2230_2249delinsGTCAA

Exon 19	Del18	c.2240_2257del
		c.2237_2255delinsT
		c.2239_2256del
		c.2236_2253del
		c.2239_2258delinsCA
		c.2237_2254del
		c.2238_2255del
		c.2237_2257delinsTCT
		c.2236_2255delinsAT
		c.2236_2256delinsATC
Exon 20	Del21	c.2237_2256delinsTT
		c.2237_2256delinsTC
		c.2235_2255delinsGGT
	Del24	c.2238_2258del
		c.2236_2256del
	T790M	c.2253_2276del
	S768I	c.2369C>T
	insG	c.2303G>T
	insASV9	c.2310_2311insGGT
	insASV11	c.2308_2309insCCCAGCGTG
Exon 21	L858R	c.2308_2311delinsCCACCGTGGAT
		c.2311_2312insGCGTGGACA
		c.2319_2320insCAC
		c.2573T>G
		c.2573_2574delinsGT
		c.2573_2574delinsGA
		c.2582T>A
	L861Q	

EGFR Total (acting as Sample Processing Control)

Specimen requirements

Sample Type	1 x 5µm FFPE tissue section
Neoplastic cells	≥10%, if less macrodissection is required
Performance	
Analytical Sensitivity	LOD ≤5% for most prevalent <i>EGFR</i> mutations
Between Laboratory	100% agreement for 10% EGFR G719S
Reproducibility	100% agreement for 10% EGFR Del15
(600 results at 3 sites)	100% agreement for 10% EGFR T790M
	100% agreement for 10% EGFR L858R
	100% agreement for 10% EGFR L861Q

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