



designed for scientists



## RSE XL Stirring bar remover

/// Data Sheet

Universal magnetic stirring bar remover for large magnetic stirring bars.

Includes mechanical coupling for removing the IKAFLON 118 float mounted stirring drive.

Length: 1.2 m

Material: PTFE, FKM (FDA-compliant)

Note: Before removing stirring bars from large vessels, please remove the I-MAG drive unit.



designed for scientists

### Technical Data

Material in contact with medium	PTFE
Depth [mm]	105
Height [mm]	1200
Weight [kg]	4

