

Cenová nabídka**Č. 22323600****Datum: 10.06.2024**

Kupující
Univerzita Palackého v Olomouci
CATRIN
Křížkovského 511/8
779 00 Olomouc

Adresát
Univerzita Palackého v Olomouci
Ústav molekulární a transl. medicíny, LF
[REDACTED]
Hněvotínská 5
779 00 Olomouc

Informace				
Vaše zákaznické č.	1256081	Kontaktní osoba	[REDACTED]	
Vaše referenční č./Datumemail_7.6.2024/		Telefon	[REDACTED]	
Platí do	10.07.2024	E-mail	[REDACTED]	
Datum dodání	11.06.2024	Koncový uživatel	Ústav molekulární a translační medi LF, Univerzita Palackého v Olomouci	
Incoterms	DAP Olomouc			
Platební podmínky	30 dnů			

Položka	Materiál/Popis	Množství	Jednotková cena CZK	Částka CZK
10	1862649	1 KS	12.700,00 (-5%)	12.700,00 -635,00
			Item Net	12.065,00
	His-Tag Capture Surface Creation Kit The HCK II contains ready-to-use solutions to create your own His-Tag Capture surface on a High Capacity Amine Sensor for use in Sierra SPR instruments. Content: 5 mL ANTA/K ₂ CO ₃ , 30 mL HEB buffer, 30 mL EDTA, pH 8.0 (100 mM), 30 mL NiSO ₄ (0.2 mM)			
20	1862614	1 KS	17.300,00 (-5%)	17.300,00 -865,00
			Item Net	16.435,00
	High Capacity Amine Sensor (3 pk.) The High Capacity Amine Sensor (HCA) contains a carboxymethylated dextran matrix covalently attached to the gold surface to extend the region where biomolecular interactions take place thereby increasing the surface capacity. It provides a hydrophilic environment and the			

BRUKER s.r.o.
Pražákova 1000/60
619 00 Brno, Czech Republic
[REDACTED]
www.bruker-sro.cz

IČO: 28297211, zapsána v OR u Krajského soudu v Brně, oddíl C,
vložka 59455
DIČ: CZ 28297211
Bankovní spojení: [REDACTED]



Položka	Materiál/Popis	Množství	Jednotková cena CZK	Částka CZK
	basis for covalent coupling of biomolecules to the surface. The negatively charged carboxyl groups on the dextran allow ligand immobilization via its functional groups, such as amino, thiol, and aldehyde groups. Alternatively, a capturing molecule may be coupled for subsequent capture. The HCA Sensor is suitable for most general applications including interaction analysis of small organic molecules, proteins, lipids, carbohydrates, and oligonucleotides. Size: 3 pk. Temperature Sensitive Product! Please store immediately at indicated temperature on arrival. Shipping at ambient temperature has been tested to be acceptable.			
30	1862617	1 KS	34.600,00 (-5%)	34.600,00 -1.730,00
			Item Net	32.870,00
	His-Tag Capture Sensor (3 pk.) The His-Tag Capture Sensor is based on a carboxymethylated dextran sensor that contains a functionalized amino-NTA surface to capture poly-His-tagged molecules through metal chelation. This sensor has the potential for efficient surface regeneration by injection of EDTA or imidazole solution which removes the metal ion and releases the His-tagged ligand. It is recommended to use the His-Tag Capture Sensors with the "His-Tag Capture Reagent Kit" (1862650) containing NiSO₄- and EDTA regeneration solutions. Size: 3 pk. Temperature Sensitive Product! Please store immediately at indicated temperature on arrival.			
40	1862634	1 KS	9.510,00 (-5%)	9.510,00 -475,50
			Item Net	9.034,50
	Amine Coupling Kit I (12-100 immob) The ACK I contains all reagents for covalent immobilization of molecules that carry a primary amine group. Sufficient for 12-100 immobilizations. Content: 750 mg EDC, 57.5 mg NHS, 20 mL Ethanolamine-HCl solution, pH 8.5. UN 2811 TOXIC SOLID, ORGANIC, N.O.S. (contains 1-(3-Dimethylaminopropyl)-3-ethylcarbodiimidehydrochloride), 6.1, III, (E), ENVIRONMENTALLY HAZARDOUS Transport category 2.			



Položka	Materiál/Popis	Množství	Jednotková cena CZK	Částka CZK
Temperature Sensitive Product! Please store immediately at indicated temperature on arrival.				
50	1862984	1 BAL	7.560,00 (-5%)	7.560,00 -378,00
			Item Net	7.182,00
	96-well (Short) Microtiter Plates 96-well Microtiter Plates (short well) compatible with all Bruker SPR instruments; working volume: 50 - 300 µl; pack of 100 plates			
60	1892529	1 BAL	16.400,00 (-5%)	16.400,00 -820,00
			Item Net	15.580,00
	96-well (Short) Microplate Sealer Self-adhesive plate sealers to cover 96-well microtiter plates. Sealers are easily pierceable and not pre-cut in order to prevent evaporation during storage. Sealers are mainly used for short well plates containing samples but can also be used for deep well plates. Compatible with polystyrene and polypropylene microtiter plates. Pack of 100 foils. Research use only.			
Netto			CZK	93.166,50
Výstupní DPH			21,00%	19.564,97
Brutto			CZK	112.731,47
Platí naše obchodní podmínky. Jsou uvedeny na adrese: https://www.bruker-sro.cz/images/obchodni-podminky.pdf				

