



SI Scientific Instruments GmbH

Römerstr. 67, 82205 Gilching, Germany

Fyzikální ústav AV ČR v.v.i.  
Na Slovance 2

182 21 Praha 8

CZECH REPUBLIC

Römerstr. 67      Telefon: 08105/77940  
82205 Gilching, Germany      Telefax: 08105/779422  
Internet: www.SI-GmbH.de      EMail: Info@SI-GmbH.de

ORDER CONFIRMATION  
NO.: AB-166-07/19  
Your order of 23.07.2019  
Mr. Ing. Urbanec Andrej  
Order No.: 0019920219  
VAT NO: CZ68378271

Datum/Date:      25.07.2019 / Mo

Dear Sirs,

we thank you for your order and confirm as follows:

Pos. Item	Menge Quant.	Bezeichnung Description	Betrag Amount
1	40	<b>MX200 HIGH PERFORMANCE VOLTAGE AMPLIFIER</b>	<b>12.400,00 €</b>
		<b>TOTAL:</b>	<b>12.400,00 €</b>
		Prices:                      Excl. local VAT. Terms of Payment:      14 days net (by bank transfer only). Delivery:                      about 4 weeks Terms of Delivery:      Ex works Gilching, Germany. Warranty:                      None, since parts or components.  Sincerely yours SI Scientific Instruments GmbH  i.A. Melanie Morstadt	

All instruments / software delivered by SI GmbH are subject to the export regulations of Germany and the USA.

Firmensitz:  
Römerstr. 67  
82205 Gilching,  
Germany  
Tel.: 08105 / 77940  
Fax: 08105 / 779422  
Email: Info@SI-GmbH.de

Büro Berlin:  
Grabenstr. 35  
12209 Berlin,  
Germany  
Tel.: 030 / 77391590  
Fax: 030 / 77391592

Geschäftsführer: Michael Dreher  
  
Handelsregister: München HRB  
69593  
USt.-IdNr. DE 128244540

Bankverbindungen:  
Raiffeisenbank Gilching  
Swift Code (BIC): GENODEF1GIL  
HypoVereinsbank München  
Swift Code (BIC):  
HYVEDEMMXXX  
HypoVereinsbank München  
Swift Code (BIC):  
HYVEDEMMXXX

€-Konto:                      62 987 (BLZ 701 693 82)  
IBAN:                              DE08 7016 9382 0000 0629 87  
48 231 900 (BLZ 700 202 70)  
€-Konto:                      IBAN:                              DE02 7002 0270 0048 2319 00  
848 227 900 (BLZ 700 202 70)  
\$-Konto:                      IBAN:                              DE05 7002 0270 0848 2279 00