

Nabídka

ÚSTAV PRO HYDRODYNAMIKU AV ČR, v.v.i.
 Pod Pařankou 30/5
 166 12 Praha 6

Informace

Číslo dokladu	4200006608	Datum dokladu	27.09.2019
Referent	Pan Holomek	Telefon	
ID účtu	1220564	DIČ	CZ67985874
Incoterms	Náklady, pojištění & přepravné	Plateb.podm.	30 dní od vystavení faktury
Podmínka expedice	Truck	Dodací lhůta	see note
Platí do	30.11.2019		

Položka	Materiál/Popis/Podmínka	Množ.	Jednotková cena	Hodnota
10	638-91112-58 TOC-LCPH FA,E200			
	Basic unit: -Measurement method:catalytic combustion at 680°C with NDIR-detection -Measurement range:TC:up to 30000 mg/L -IC :up to 35000 mg/L -Detection limit: 4 µg/L -Determination of the parameter TC, TIC, TOC (with difference method) and NPOC -Blank-Check-Function to determine the blank value -Normal sense catalyst with an injection volume of maximal 150µL -High sense catalyst with an injection volume of maximal 2000 µL -Simple change of catalyst -Automatic injection with the well known syringe/slider technique (without septa) -ISP-module (integrated sample pre-treatment): -Automatic dilution from factor 2 to 50 in the syringe -Acidification and sparging for NPOC-method in the closed environment			

Str. 1 / 5

Nabídka

4200006608

Položka Materiál/Popis/Podmínka	Množ.	Jednotková cena	Hodnota
<p>(free from environmental contamination)</p> <ul style="list-style-type: none"> -Multi-point calibration from one Stock solution -Optimised homogenisation of suspended particle due to high sampling volume -Automatic change of measurement parameter if the sample exceed the measurement range -Automatic calibration -EN 1484 conform <p>Options:</p> <ul style="list-style-type: none"> -TNM-L, upgrade to measure Total Nitrogen (TNb) -POC option, upgrade to addition method (TOC=NPOC+POC) -SSM-5000A, upgrade to measure TOC in solids -High Salt Kit to measure salted samples -Small sample volume kit for small sample volumes -Manual injection kit, <p>Operation - PC-controlled</p> <p>TOC-Control L for TOC-L CPH to control the system completely via USB connection:</p> <p>User Friendly Interface</p> <ul style="list-style-type: none"> -Icons and functions are shown on large buttons -Lists files being used by types, sorting by file name or creation date -Samples can be inserted in the sample table by dragging method files -Select operations from the right-click menu <p>Useful functions facilitate your analysis operations</p> <ul style="list-style-type: none"> -Edit mode allows adding of samples, printing or recalculation of measured samples during measurement -Create schedule files for routine analysis -Recalculation of several samples with another calibration curve -User defined data report -Various Export/Import functions -Adaptation to a LIMS-system <p>Quality assurance functions and GLP/GMP conform</p> <ul style="list-style-type: none"> -Control samples, control charts -Data management including raw data -Individual User Administration -Change of Log in during measurement -Audit trail <p>Supports 21 CFR Part 11 -</p> <p>Together with Shimadzu's LabSolutions database, the TOC data can be acquired, archived and secured. The electronic signature can be used to release the secured data (LabSolutions).</p>			

Nabídka
4200006608

Položka	Materiál/Popis/Podmínka	Množ.	Jednotková cena	Hodnota
		1,00 KS		
20	638-93199-58 ASI-L FA 24 E			
	ASI-L 24ml Autosampler base unit for TOC-L series	1,00 KS		
30	638-92325-41 ACCESSORIES, ASI-L 24			
	Accessory set for ASI-L 24ml 93 position sample rack, sampling needle, communication cable	1,00 KS		
40	638-41462 VIAL SET 24ML. 100 PCS	2,00 KS		
50	038-00165-61 SEPTUM FOR VIAL 24 ML. 100 PCS	2,00 KS		
60	638-20074-01 CAP. FOR VIAL 24ML. 100PCS	2,00 KS		
70	638-67099-41 STIRRER ASSY,24A L			
	Internal magnetic stirrer option for ASI-L 24ml (sample positions 1 to 85 of 93)	1,00 KS		
80	046-00617-02 STIRRER BAR,SA-10			

Str. 3 / 5

Nabídka
4200006608

Položka	Materiál/Popis/Podmínka	Množ.	Jednotková cena	Hodnota
		20,00		
90	638-77183-41 EXTERNAL SPARGE ASSY,ASI-L			
	Allows sparging inside ASI-L sample vials	1,00		
100	638-91108-58 TNM ASSY,L			
	Upgrade option to measure Total Nitrogen Measurement method: <(>,<)>catalytic combustion at 720°C with Chemiluminescence detection EN 12260 conform TNb and TOC can be measured simultaneously Measurement range: <(>,<)>up to 10000 mg/l Detection limit: <(>,<)><(>,<)>5µg/l (for TOC-L CPH, CSH) <(>,<)><(>,<)><(>,<)><(>,<)><(>,<)>20µg/l (for TOC-L CPN, CSN)	1,00		
110	638-56255-41 CARRIER GAS CONTROL VALVE SET			
	Internal carrier gas pressure regulator (required inlet pressure 3 - 6bar)	1,00		
120	630-00992 HALOGEN REMOVER CARTRIDGE			
		1,00		
130	638-41323 TC COMBUSTION TUBE WITH CASE			
		1,00		
140	630-00999 CO2 ABSORBER			
		1,00		

Str. 4 / 5

Nabídka

4200006608

Položka	Materiál/Popis/Podmínka	Množ.	Jednotková cena	Hodnota
150	630-00105-01 PLATINUM NET 2 PCS.TOC-500	1,00 KS		
160	INSTALLATION- INSTALLATION TOC TOC	2,00 KS		
170	TRAINING-TOC TRAINING TOC	1,00 KS		
180	638-60116 CATALYST PT.ST TYPE	1,00 KS		
190	980-23765 HP ProDesk 600 G4 MT i5-8500/8GB/1TB/ DV	1,00 KS		
200	980-23770 HP E233 23" IPS 192x1080/250/1000:1/ VGA	1,00 KS		
Cena bez DPH				1.437.913,60
Výstupní DPH 21,00 %				301.961,86
Cena s DPH				CZK 1.739.875,46

Cena zahrnuje dopravu na místo instalace, instalaci přístroje a odborné školení uživatelů.

Záruční doba je 24 měsíců od instalace, maximálně 26 měsíců od data dodání.

Platební podmínky - fakturace po instalaci, splatnost 30 dní není-li dohodnuto jinak v kupní smlouvě.

V ceně je účtován poplatek za likvidaci elektroodpadu dle zákona 185/2001Sb. vedenému v kolektivním systému RETELA pod registračním číslem 2005/10/10/439

SHIMADZU Handels GmbH organizační složka