

Fyzikální ústav Akademie věd ČR

Na Slovance 1999/2
182 21 Praha 8

Quote No.: 49229-276840
Customer No.: 061001109

Date: 2.11.2021

Your Struers Rep.: [REDACTED]

Re.: LaboPol-30 + UI

[REDACTED]

Related to your request, we would like first to thank you for your interest in our products and we hereby send you a price quotation for the required equipment:

- LaboPol-30 with a manual control panel UI and magnetic disc.

Should you have any further questions related to that quotation, please, do not hesitate to contact me. Thank you and we are looking for our further cooperation.

Kind regards,






STRUERS GmbH, organizational unit

[REDACTED]
Sales Engineer
[REDACTED]

Attachments: (attachments are an integral part of that quotation)

1) Information brochure to LaboSystem.



Cat. No.	Product description	Qty.	Unit Price	Total CZK
06336127	LaboPol-30  Robust manual grinding and polishing table-top 250 or 300 mm dia. disc solution for fast, reliable specimen preparation 24/7. With a flexible modular design that is simple to scale-up with a specimen mover if your requirements increase. LaboSystem main features: - Fast - Shortest possible turnaround - Reliable - Tested and proven to withstand hard work - Modular system - Different units can be combined and the system can be upgraded later - Removable bowl liner - Splash ring for every purpose Supplied with splash guard for manual preparation with a 300 mm dia. disc and bowl liner. LaboUI control panel (06206901) or specimen mover LaboForce-50 (06356127)/LaboForce-100 (06366127) and discs are ordered separately. 1/3 x 200-240 V / 50-60 Hz	1	128.610,00	128.610,00
06206913	Lid for LaboPol-30/-60 For protection of LaboPol-30/-60 disc during non-use 	1	2.780,00	2.780,00
06086403	MD-Disc with cone MD-Disc for use on Tegramin-30 and LaboPol-30/-60. With cone for easy replacement 300 mm (12") dia. 	1	23.580,00	23.580,00
06206901	LaboUI  Robust control panel for manual preparation for fast preparation 24/7. LaboUI main features: - For manual preparation - Elevated placing For use with LaboPol-20/-30/-60 and LaboDoser-10 dosing unit.	1	5.590,00	5.590,00
49900049	MD-Gekko Adapter for use with SiC Foil or self adhesive consumables on the MD-System, for easy removal. For magnetic fixation on MD-Disc. 300 mm dia. 2 pcs. 	1	2.849,00	2.849,00
STARTUP1	Start-up Package - Product group 1	1	6.486,90	6.486,90

Cat. No.	Product description	Qty.	Unit Price	Total CZK
	Struers installs the equipment and prepares it for start-up. We provide basic hands-on training in the main function of your equipment. Your staff is trained in daily, weekly & monthly maintenance. A dedicated application training is done to make sure the correct method is used. Travel and accommodation included. 2 hours			

Sub Total: 169.895,90

Total excl tax 169.895,90
 Tax (21,00 %) 35.678,14
Total incl tax 205.574,04

Terms & Conditions

Valid until:	30.12.2021
Shipment terms:	CIP
Payment terms:	30 days
Delivery time:	6-8 weeks
Warranty in years:	2

General Terms & Conditions: www.struers.com/Terms-and-Conditions-of-sale

If you want to receive Struers information about training, free webinars and new products please sign up here:
www.struers.com/en/Academy/Newsletter



When teaming up with Struers

You optimize processes and results

Improve methods and workflows – with specialist consulting and troubleshooting

Boost your knowledge – webinars, articles and application notes keep you up to date

Increase your skills – training courses and workshops give you hands-on understanding of the best processes and techniques for your lab

You protect your investment

Improve efficiency – create optimized solutions by combining our products for highest efficiency and productivity

Minimize downtime – with high quality equipment, maintenance services and support hotlines

Enjoy long-term performance – we supply spare parts for at least 10 years after equipment has been discontinued

You maximize uptime with fast delivery

Standard shipment – of consumables and spare parts within 48 hours

No waiting time – your orders arrive complete and shipped on the promised date

Keep your stock costs down – 98% of our consumables and spare parts are kept in stock

You guard your working environment

Experience valid and reproducible results – our ISO 9001-certified LEAN production and strict quality control ensure your equipment and consumables will perform as promised

Be safe – we follow product classifications ensuring the safest possible working environment

Minimize environmental impact – with ISO 14001 certification and low-carbon shipping