


Swing-out rotor 11118

with 1 set (qty 2) rectangular buckets 13218, caps 17108 and plate holder 17978

- Max. plate height without / with cap: 93 / 80 mm
- Max. speed (6-16S / 6-16KS): 5,700 rpm
- Max. RCF (6-16S / 6-16KS): 5,485 x g
- Max. load per carrier: 520 g
- Biosafe with cap 17108 
- Tmin at maximum speed: <4°C
- Angle: 90°



Swing-out rotor 11660 / 11650*

Rectangular bucket 13421 (1 set = 4 pcs) with plate holder 17979

- Max. plate height: 68 mm
- Max. speed (6-16S / 6-16KS): 4,200 / 5,100 rpm
- Max. RCF (6-16S / 6-16KS): 3,096 / 4,565 x g
- Max. load per carrier: 500 g
- Tmin at maximum speed: <4°C
- Angle: 90°



*Max. speed with rotor 11650 in 6-16KS: 4,700 rpm, max. RCF: 3,877 x g

Plate



Description

MTP (without lid)	MTP (with lid)	DWP	PCR plate 96-fold	Filter plate
-------------------	----------------	-----	-------------------	--------------

W x D x H [mm]

128 x 86 x 15	128 x 86 x 17.5	128 x 86 x 44.5	128 x 86 x 20	128 x 86 x 60
---------------	-----------------	-----------------	---------------	---------------

Adapter

for rectangular bucket 13218



Bore diameter x length [mm]

		–		
--	--	---	--	--

Plates per rotor (without / with cap)

12 / 10	10 / 8	4 / 2	8 / 8	4 / 4
---------	--------	-------	-------	-------

Centrifuging radius [mm]

		151		
--	--	-----	--	--

Item no.

		17978		
--	--	-------	--	--

Plate



Description


MTP (without lid)	MTP (with lid)	DWP	PCR plate 96-fold	Filter plate
-------------------	----------------	-----	-------------------	--------------

W x D x H [mm]

128 x 86 x 15	128 x 86 x 17.5	128 x 86 x 44.5	128 x 86 x 20	128 x 86 x 60
---------------	-----------------	-----------------	---------------	---------------

Adapter

for rectangular bucket 13421



Bore diameter x length [mm]

		–		
--	--	---	--	--

Plates per rotor

16	12	4	12	4
----	----	---	----	---

Centrifuging radius [mm]

		157		
--	--	-----	--	--

Item no.

		17979		
--	--	-------	--	--