

NABÍZENÁ KONFIGURACE

I. EU SIGNA Voyager AIR IQ Edition 29.1

MN.	Popis
1	SIGNA™ VOYAGER 1.5T 33CH MR System - u
1	SIGNA™ Voyager 33 to 49 Channel Upgrade
1	SIGNA™WORKS AIR™ EDITION MR29.1 - c
1	COMPUTING PLATFORM AND DICOM CONFORMANCE - T5820 GOC - c
1	Gen 7 DL Performance ICN - c
1	AIR™ Recon DL and AIR x™ Package - Voyager
1	Comfort Plus Fixed Patient Table Collector - AIR™ Edition
1	SIGNA Voyager 1.5T Magnet Collector
1	HYPERWORKS
1	DISCO Trio
1	DIFFUSION
1	Cardio Plus Package
1	oZTEo for Silent Suite
1	3D ASL
1	3D PROMO
1	Silent Suite
1	IDEAL IQ
1	Flow Analysis 4.0
1	CardioMaps
1	1.5T 16-channel AIR Anterior Array (AA)
1	1.5T AIR™ Multi-Purpose Coil Large with Positioners
1	1.5T 16-channel T/R Knee Array
1	1.5T 16-Channel T/R Hand-Wrist Coil

II. EU AW Server 3.2 for MR modality

MN.	Popis
3	Volume Illumination
3	3D Suite
1	ReadyView NRR
1	Body View
1	Brain View
1	Quantib Brain Package