



Objednávka č.: 06/2017
V Praze dne 26.7.2017

KONEKO marketing, spol. s r.o.
pí. Jiřina Neužilová
Na Cikánce 2
153 00 Praha 5 - Radotín

Věc:
Objednávka

Vážená paní Neužilová,
na základě Vaší nabídky u Vás objednávám chladicí termostat Julabo F25-ME v celkové ceně 129.833,-Kč včetně DPH.
Nedílnou součástí objednávky č. 06/2017 je Vaše nabídka.

Prosím, o potvrzení objednávky, konečné celkové ceny a termínu dodání.

Děkuji a těším se na spolupráci.

S pozdravem


Eva Krejčíková
grantová sekce

Příloha:
Nabídka na termostat F25-ME

ÚSTAV CHEMICKÝCH PROCESŮ AV ČR, v.v.i.
165 02 Praha 6 - Suchbát, Rozvojová 135
IČO: 67985858 DIČ: CZ67985858
- 3 -

Objednatel:
Ústav chemických procesů AV ČR, v.v.i.
Rozvojová 135
165 02 Praha - Suchbát
IČ: 67985858
DIČ: CZ67985858

PRÍLOHA:

KONEKO marketing, spol. s r. o.

Fakturační adresa: Sojovická 2, 197 00 Praha 9 – Kbely

Obchodní a poštovní adresa: Na Cikánci 2, 153 00 Praha 5 – Radotín

Tel./fax: 257 911 117 | mobil: 603 357 809

e-mail: koneko@laboratore.cz



www.laboratore.cz

e-shop: www.obchod.laboratore.cz

Petr Stanovsky Department of Multiphase Reactors
Institute of Chemical Process Fundamentals
Czech Academy of Sciences
Rozvojova 135/1
16502 Prague 6 - Suchbátov

Tel. +420 220 390 162

E-mail: stanovsky@icpf.cas.cz

<http://www.icpf.cas.cz>

12.4.2017

Věc: Nabídka na termostát

chladicí termostát Julabo F25-ME



Kód zboží:

F25-ME

Výrobce:

Julabo

Cena bez DPH:

107 300,00 Kč

Cena s DPH:

129 833,00 Kč

DPH:

21,00 %

Doprava : 500,- Kč + DPH

tepl. rozsah	tepl. stabilita	volný pracovní
[°C]	[°C]	prostor vany
-28 - +200	± 0,01	120 x 140 x 140

- kromě typu ED mají rozhraní RS 232
- výkonné chladicí lázně se zabudovaným termoregulátorem ED
- možnost externí cirkulace média mimo lázeň
- maximální tlak pumpy 0,34 bar, průtok 15 litrů za minutu
- objem vany 4,5 l

otvor lázně / hloubka	12 x 14 / 15 cm
výkon tlakové pumpy	11 - 16 l/min; 0,45 bar
chladicí spirála	chladicí term.

zapsaná v OR u Městského soudu v Praze

odd. C, vložka 6314

IČO: 44792964

DIČ: CZ44792964

ČSOB hlavní pobočka Praha

č. účtu: 579311683/0300

KONEKO marketing, spol. s r. o.

Fakturační adresa: Sojovická 2, 197 00 Praha 9 – Kbely

Obchodní a poštovní adresa: Na Cikánci 2, 153 00 Praha 5 – Radotín

Tel./fax: 257 911 117 | mobil: 603 357 809

e-mail: koneko@laboratore.cz



www.laboratore.cz

e-shop: www.obchod.laboratore.cz

Ceny jsou platné do konce června 2017

Dodací lhůta : cc 3-4 týdny

Placení fakturou

Další dotazy Vám rádi zodpovíme.

Jiřina Neužilová

KONEKO marketing, spol. s r.o.

Innovative Temperature Technology

JULABO



JULABO F25-ME Refrigerated and Heating Circulator

Refrigerated and Heating Circulator for internal and external temperature applications

JULABO Refrigerated & Heating Circulators for heating and cooling are suitable for external temperature tasks and for temperature applications directly in the bath. A model from three different series' with different capacities can be selected. Powerful circulating pump systems and high heating and cooling capacities guarantee short heat-up and cool-down times. The inevitable dust accumulation that occurs in refrigeration systems is solved by means of an easily removed venting grill, which allows for easy cleaning. Depending on the model, the units are equipped with handles or castors for an easy transportation. All models have a drain tap on the front allowing the fluid to be drained easily. The compact assembly of each unit enables a space to be saved in the laboratory.



Available downloads:

- + [EasyTemp Software](#)
- + [Circulators/Baths catalog \(PDF, English\)](#)
- + [Interactive product demo](#)
- + [Easy BlackBox remote diagnosis software](#)

Features

- o VFD COMFORT DISPLAY, resolution 0.1 °C
- o Keypad for setpoints, warning/safety values and menu functions
- o PID3 cascade temperature control
- o ATC3 3-Point-Calibration
- o Pt100 External sensor connection for measurement and control
- o PPC (Pump Pressure Control), electronically adjustable pump capacity
- o Early warning system for low liquid level (DBGM 203 06 059.8)
- o Adjustable high temperature cut-out, visible via display
- o RS232 interface for online communication
- o Integrated programmer for 10 program steps
- o Active Cooling Control

Model	F25-ME			
JULABO Order No.	9160625			
Working temp. range	°C	-28 ... 200		
Temperature controller		PID, cascade		
Temperature stability	± °C	0.01		
Display		VFD		
Display resolution	°C	0.1		
Integrated Programmer		yes		
Heater capacity	W	2000		
Cooling capacity	°C	20	0	-20
	W	260	200	60
Refrigerant		R134a		

Pump

Pressure pump capacity max.	mbar	450
	l/min	11-16
Pump connections		M10 x 1

Electrical connections

Digital interfaces	RS232	RS485	PROFIBUS
	X	-	-
Analog interfaces	Programmer input		-
	Temperature recorder output		-
	Stand-by		-

	Alarm	-
	External Pt100 sensor	X

General information

Bath opening (W x L)	cm	12 x 14
Bath depth	cm	14
Filling volume	Liters	4.5
Dimensions (W x L x H)	cm	23 x 42 x 61
Weight	kg	31
Ambient temperature	°C	5...40
Mains power (±10%)	V/Hz	230V / 50Hz or 115V / 60Hz or 100V / 50-60Hz or 230V / 60Hz
Power consumption	W	n/a
Classification according to DIN 12876-1	Class	3 (FL)

Miscellaneous

Removable venting grid for convenient cleaning to maintain cooling performance.

Data determined at a mains voltage of 230V / 50Hz and ambient temperature of 20 °C.
Technical changes without prior notification reserved. Pictures may differ from original.

JULABO Labortechnik GmbH
D-77960 Seelbach / Germany
Telephone +49 7823/51-0 | Telefax +49 7823/2491
Email info@julabo.de